

The right temperature worldwide

# LAUDA



## Chapter Microcool

### LAUDA – the big one

Thermostats, Circulation chillers, Water baths

Overall Brochure 2016/2017

# LAUDA Microcool

Circulation chillers for reliable continuous operation in the lab and in research from -10 up to 40 °C



Excellent price-performance ratio, compact design and simple to use



## Application examples

- Cooling of rotary evaporators
- Cooling of distillation systems
- Supply of cooling traps
- Cooling of analytical devices

**LAUDA Microcool** has been designed as a circulation chiller line with five compact models and cooling capacities from 0.25 to 1.2 kW. The user interface with large LED display and the membrane keyboard make the devices easy to use. An RS-232 interface and alarm contact are integrated as standard. What is unusual in this price category

of circulation chillers is the high-quality block pump with magnetic coupling. The magnetic coupling of pump and electric motor exclude sealing problems on the pump shaft. LAUDA Microcool circulation chillers are used whenever heat needs to be dissipated reliably and quickly, e.g. in laboratories for rotary evaporators, distillation systems or analytical devices.



# Your advantages at a glance

+	The Microcool advantages	Your benefits
	<ul style="list-style-type: none"> <li>• Five device types in four housing sizes</li> <li>• Cooling capacities from 250 W up to 1200 W</li> </ul>	<ul style="list-style-type: none"> <li>• Clear device portfolio for simple selection</li> <li>• Covers the majority of basic lab uses</li> </ul>
	<ul style="list-style-type: none"> <li>• User interface with large LED display and membrane keyboard</li> <li>• Autostart timer and auto-shutdown function</li> <li>• Illuminated window for checking heat transfer liquid level</li> </ul>	<ul style="list-style-type: none"> <li>• Simple and intuitive use</li> <li>• Timer-based activation and deactivation of the circulation chillers</li> <li>• Quick optical detection of the filling level</li> </ul>
	<ul style="list-style-type: none"> <li>• Block pump with magnetic coupling of pump and electric motor</li> <li>• Integrated adjustable bypass and pressure gauge at MC 600, MC 1200 and MC 1200 W</li> <li>• Integrated overflow connection</li> </ul>	<ul style="list-style-type: none"> <li>• Prevents sealing problems at the pump shaft</li> <li>• Integrated pump pressure adjustment for connected delicate glassware</li> <li>• Controlled filling of the devices</li> </ul>
	<ul style="list-style-type: none"> <li>• RS-232 interface and alarm contact standard</li> </ul>	<ul style="list-style-type: none"> <li>• System integration into processes without additional costs</li> </ul>
	<ul style="list-style-type: none"> <li>• Compact design and low space requirements</li> <li>• Integrated filling funnel on top of the device</li> <li>• Easily removable front grid</li> </ul>	<ul style="list-style-type: none"> <li>• Saves valuable laboratory space</li> <li>• Simple and safe filling</li> <li>• Easy-to-clean condenser</li> </ul>

# LAUDA Microcool

## Microcool Circulation chiller with cooling capacity up to 1200 Watt

The compact MC 250 and MC 350 makes them ideal for being positioned on the benchtop. The circulation chillers are equipped with a magnetic coupling pump. This supplies a pump pressure of 0.35 bar and a maximum pump flow of 16 L/min.


The 600 and 1200 Watt cooling capacity models are floor standing instruments designed to fit underneath the lab bench. They are equipped with a pressure gauge to display the pressure and casters which can be controlled and locked. Pump pressure can be adjusted via the integrated bypass. At 1200 Watt, the most powerful device is also available in a water-cooled version as the MC 1200 W.



Circulation chiller MC 250



### Pump characteristics Heat transfer liquid: Water



- ① MC 250, MC 350
- ② MC 600  
MC 1200  
MC 1200 W

### Temperature range

-10...40 °C

### Included as standard

RS 232 interface · alarm contact

### Included accessories (except of MC 250, MC 350)

Nipples (3/4") · screw caps



All technical data on page 106  
and following

Other power supply variants on page 114



465 mm



500 mm



595 mm



650 mm



650 mm

Technical features		MC 250	MC 350	MC 600	MC 1200	MC 1200 W
Working temperature range*	°C	-10...40	-10...40	-10...40	-10...40	-10...40
Temperature stability	±K	0.5	0.5	0.5	0.5	0.5
Cooling output at 20 °C	kW	0.25	0.35	0.6	1.2	1.2
Pump pressure max.	bar	0.35	0.35	1.3	1.3	1.3
Pump flow max.	L/min	16	16	35	35	35
Cat. No. 230 V; 50 Hz		LWM 118	LWM 119	LWM 120	LWM 121	LWM 122

\* Working temperature range is equal to ACC range

## Microcool accessories (excerpt)

### EPDM tubing

Cat. No.	Description	d <sub>i</sub> (mm)	d <sub>e</sub> (mm)	Temperature range °C	Pressure range max. bar
RKJ 111	Polymer tubing	9	11	10...120	1
RKJ 112	Polymer tubing	12	14	10...120	1
LZS 021	Insulated	12	21	-35...90	-
RKJ 031	Reinforced fibres	13 (1/2")	19	-40...100	20
RKJ 032	Reinforced fibres	19 (3/4")	27	-40...100	20
RKJ 009	Tube insulation for 1/2"	23	33	-50...105	-
RKJ 013	Tube insulation for 3/4"	29	39.5	-50...105	-

d<sub>i</sub> = internal diameter ; d<sub>e</sub> = external diameter

### Adapter G 3/4"

Cat. No.	Designation	Description
LWZ 016	Nipple	3/4" Screw cap, 1/2" nipple
LWZ 040	Nipple	3/4" Screw cap, 10 mm nipple

### Stainless steel hose clamps

To secure hoses

Cat. No.	Description
EZS 012	Hose clamp for external diameter 10-16 mm, 1/2"
EZS 013	Hose clamp for external diameter 12-22 mm, 1/2"
EZS 015	Hose clamp for external diameter 20-32 mm, 3/4"

### Heat transfer liquids

Cat. No.	Description	Temperature range °C
LZB 120	Aqua 90, 5 L	5...90
LZB 220	Aqua 90, 10 L	5...90
LZB 320	Aqua 90, 20 L	5...90
LZB 109	Kryo 30, 5 L	-30...90
LZB 209	Kryo 30, 10 L	-30...90
LZB 309	Kryo 30, 20 L	-30...90



RKJ 031



LWZ 016



EZS 012



LZB 209



Order the detailed LAUDA accessories brochure and the heat transfer liquids brochure free of charge. These and additional product information can also be found at [www.lauda.de](http://www.lauda.de)

## Temperature solutions:

Thermostats · Circulation chillers · Water baths

Process cooling systems · Heat transfer systems · Secondary circuit systems



**LAUDA DR. R. WOBSE  
GMBH & CO. KG**  
Headquarters  
Pfarrstraße 41/43  
97922 Lauda-Königshofen  
Germany  
Phone: +49 (0)9343 503-0  
E-mail: [info@lauda.de](mailto:info@lauda.de)



**LAUDA-Noah, LP**  
308 Digital Drive  
Morgan Hill, CA 95037  
USA  
Tel.: +1 360 993 1395  
E-mail: [info@lauda-noah.com](mailto:info@lauda-noah.com)



**LAUDA Technology Ltd.**  
4200 Waterside  
Solihull Parkway  
Birmingham Business Park  
B37 7YN Birmingham  
Great Britain  
Phone: +44 121 717 4789  
E-mail: [info@lauda-technology.co.uk](mailto:info@lauda-technology.co.uk)



**LAUDA China Co. Ltd.**  
**Shanghai**  
2nd floor, Building 6  
No. 201 MinYi Road  
Songjiang District  
201612 Shanghai  
China  
Phone: +86 21 64401098  
E-mail: [info@lauda.cn](mailto:info@lauda.cn)



**LAUDA-Brinkmann, LP**  
1819 Underwood Boulevard  
08075 Delran, NJ  
USA  
Phone: +1 856 7647300  
E-mail: [info@lauda-brinkmann.com](mailto:info@lauda-brinkmann.com)



**LAUDA América Latina  
Tecnologia Ltda.**  
Av. Paulista, 726 – 17º andar – Cj. 1707  
01310-910 – São Paulo – SP  
Brazil  
Phone: +55 11 3192-3904  
E-mail: [info@lauda.net.br](mailto:info@lauda.net.br)



**LAUDA France S.A.R.L.**  
Parc Technologique de Paris Nord II  
Bâtiment G  
69, rue de la Belle Etoile  
BP 81050 Roissy en France  
95933 Roissy Charles de Gaulle Cedex  
France  
Phone: +33 1 48638009  
E-mail: [info@lauda.fr](mailto:info@lauda.fr)



**Office Beijing**  
15/F, Office Building A,  
Parkview Green,  
9 Dongdaqiao Road,  
Chaoyang District  
100020 Beijing  
China  
Phone: +86 10 57306210  
E-mail: [info@lauda.cn](mailto:info@lauda.cn)



**LAUDA-Brinkmann, LP**  
308 Digital Drive  
Morgan Hill, CA 95037  
USA  
Phone: +1 856 7647300  
E-mail: [info@lauda-brinkmann.com](mailto:info@lauda-brinkmann.com)



**LAUDA Ultracool S.L.**  
C/ Colom, 606  
08228 Terrassa (Barcelona)  
Spain  
Phone: +34 93 7854866  
E-mail: [info@lauda-ultracool.com](mailto:info@lauda-ultracool.com)



**LAUDA Italia S.r.l.**  
Strada 6 – Palazzo A – Scala 13  
20090 Assago Milanofiori (MI)  
Italy  
Phone: +39 02 9079194  
E-mail: [info@lauda-italia.it](mailto:info@lauda-italia.it)



**LAUDA Singapore Pte. Ltd.**  
25 International Business Park  
#04-103M German Centre  
Singapore 609916  
Phone: +65 6563 0241  
E-mail: [info@lauda.sg](mailto:info@lauda.sg)



**LAUDA-Noah, LP**  
2501 SE Columbia Way, Suite 140  
Vancouver, WA 98661  
USA  
Tel.: +1 360 993 1395  
E-Mail: [info@lauda-noah.com](mailto:info@lauda-noah.com)



**LAUDA IBÉRICA SOLUCIONES  
TÉCNICAS, S.L.**  
C/ Colom, 606  
08228 Terrassa (Barcelona)  
Spain  
Phone: +34 93 7854866  
E-mail: [info@lauda-iberica.es](mailto:info@lauda-iberica.es)



**OOO „LAUDA Wostok“**  
Malaja Pirogovskaja Str. 5  
119435 Moscow  
Russia  
Phone: +7 495 9376562  
E-mail: [info@lauda.ru](mailto:info@lauda.ru)

**LAUDA DR. R. WOBSE GMBH & CO. KG**  
Pfarrstraße 41/43 · 97922 Lauda-Königshofen · Germany  
Phone: +49 (0)9343 503-0 · Fax: +49 (0)9343 503-222  
E-mail: [info@lauda.de](mailto:info@lauda.de) · Internet: [www.lauda.de](http://www.lauda.de)