The right temperature worldwide

LAUDA



LAUDA Measuring instrumentation:

Viscosity measuring system PVS

1998



Existing PVS/VRM complemented by fully automatic sampler, becomes the VAS 1. The logical step to full automation of viscosity measuring.

1999



The improved Drop Volume Tensiometer TVT 2, equipped with the latest microprocessor technology and easy to use Windows program, was presented this year.

2000



Introduction of the new generation of Ring/Plate Tensiometer TE 2 at ACHEMA 2000.

2004



The LAUDA viscosity measuring system DVS 1 has been optimised for the parallel determination of viscosity and density. The new measuring method and uncomplicated handling enable shorter measuring times than is the case with conventional viscometers.

2005



The LAUDA vicosity measuring system PVS participates in the advantages of the LAUDA Proline thermostats. Measurements of low temperatures down to -60 °C are available.

Product advantages

Modularity, the basic principle 4
Modules at a glance 6
Autosampler VAS 1 7
Software 8

Page

9

Applications

Compatibility with laboratory environment 9
Sovent viscosity of plastics 10
Testing technical lubricants 12
Integrated solvent recycling 13
Determining enzyme activity 14

Technical Data & Accessories 15

Modules15Automatic cleaning17Viscometer for offline cleaning18Viscometer for online cleaning19Clear-view thermostats20Through-flow coolers21Equipment combinations22

Glossary

23

Modularity, the basic principle

Increasingly stringent quality specifications demand more and more accurate control of raw materials and intermediate products. Kinematic viscosity is an important characteristic of liquids with nearly Newtonian (i.e. ideal) flow behaviour. And capillary viscometry is the most accurate method for measuring it. Measurement routines should be efficient, rapid, absolutely reliable and be reproducible without limitation.

Viscosity with the flow:

With LAUDA's unique modular concept it is possible to set up system configurations which represent the optimal combination of all functionalities required for a particular application. These range from 1-place measuring systems up to 8-place systems with automatic cleaning and the 4-place system with auto-sampler. A very high degree of automation is achieved for repeated measurement routines. The numerous operations which are often still performed manually today are thus reduced to a minimum.

Using these individual configurations, viscosities and characteristics derived from them can be evaluated for a great variety of materials:

Plastics: measuring relative, reduced and intrinsic viscosity as a measure of mean molecular weight and thus of polymer length which defines quality.

Lubricants, oil and fuel: measuring viscosity and its variation with temperature as well as viscosity indexes of mineral oil products, additives and their mixtures according to ASTM and ISO standards.

Enzyme activity: determining the reaction activity of certain enzymes from their effect on the time change in the flow characteristic of dissolved biological molecules.

Cellulose: determining the chain length of basic materials for technical papers and textiles, and the change in polymer length over various processing stages and due to wear.






Modularity, the basic principle

The LAUDA viscosity measuring system is built up from independent and self-contained functional units. These are linked through a central control unit to a conventional PC which controls the entire measurement sequence and evaluates the measured data. The decentralized structure enables all compo-

nents to operate independently of each other. Through this independence of the modules it becomes possible to assign different measurement routines to individual places, so that routines can run independently on the different places without any mutual interaction.

Module adjustment to application



The modular structure offers application-related solutions with extensive extension options with regard to the degree of automation and ease of evaluation. The affordable single-place measuring system in its minimally-configured form already contains the basic components. This can be further extended through to the four-place fully-automatic systems with sample charging and automatic cleaning and integrated solvent recycling and inclusion of the sample preparation. The pressing objective of this is to make the measurements and their preparation independent of the relevant user, and to simultaneously free the user from routine activities, e.g. dealing with hazardous solvents.





Flexibility


- Individual combination of components to suit the current task
- No ballast through unnecessary functions
- Readily extended at any time
- Always the latest version by integrating newly developed modules

Functionality

- ◆ Decentralized arrangement
- Independent intelligent components through single chip processors
- Simple and reliable operation
- Long life, all functional components have maximum resistance to chemicals and heat

Modules at a glance

The modular structure also permits economic system arrangements which are optimally matched to actual requirements. The configuration can thus be adapted to a larger sample throughput, to new tasks, or to the integration of newly developed modules.



Control unit PVS 1

is the central module of the system and at the same time forms the link between the PC and individual components. The control unit provides a total of four slots which, depending on configuration, can be fitted with modules for one (ME 1) or two (ME 2) measuring stands, dosing systems (BE), or the autosampler (VAS/E) and magnetic stirrer control (MRE).



Measuring stand S 5

can carry different standard capillary viscometers, for example the types Ubbelohde (see illustration) or Cannon-Fenske Routine. The time for the sample to flow through the viscometer capillary is measured to the nearest millisecond, using a novel infrared sensor controlled by a single-chip processor. The sturdy micro pump for transferring the sample up to the bulb, together with the chemical-resistant valves in the stand head, ensure very compact construction and reliable long-term operation.



Autosampler VAS 1

U+ VAS 1 complements a 4-place measuring system to provide maximum automation convenience with an extremely high sample throughput. Up to 63 samples, depending on reservoir size, can be processed in one setting. Also hot polymer solutions or oils can be handled in heated rack with heated dosing syringe.



Cleaning module VRM








provides fully automatic cleaning and drying of the viscometers. Either one or two viscometers can be connected and two different cleaning liquids can be selected separately. Even very hot samples up to 180 °C (VRM 3D/HT) and high viscous samples (VRM 3D/HV) can be handled reliably. Use of high-grade materials ensures absolute chemical resistance.




Dosing system



dosimat for determining limiting viscosity through different concentration steps, in conjunction with a dilution viscometer, a magnetic stirrer operating from the burette module (BE) it can be controlled via measuring software.



Thermostat



Precise measurement of viscosity demands that the test temperature is kept constant and uniform throughout the bath. LAUDA clear-view thermostats, or LAUDA Ecoline Staredition thermostats in conjunction with a transparent bath are important elements permitting unrestricted observation of the capillary viscometers.



The VAS 1 autosampler is the ideal modular complement when large sample throughputs and regularly repeated measurement routines are involved. Monotonous and time consuming operations are automated, the operators can concentrate on more important tasks. The effort required per measurement, from filling the

viscometer up to its final cleaning, is greatly reduced, permitting efficient loading of the measuring system. Typically up to 160 samples can be evaluated in 10 hours, and correspondingly more with shift operation. The operator no longer comes into contact with irritating, corrosive, poisonous or hot substances.

Better reproducibility and precision



The autosampler fills up to four viscometers in parallel with the same or different samples. A syringe draws up the sample from closed reservoirs which are mounted on a sample stand with optional thermostating capability. The syringe is automatically cleaned to prevent carry-over of a previous sample or cleaning solvent into the sample on test, a common danger with the more usual fixed tubing connections. The result is an appreciable improvement in reproducibility and accuracy. Before dosing the next sample, the VRM modules ensure thorough cleaning and drying of viscometers and dosing syringe.


The sequence of the individual tests and the sample assignment are determined by the PC or can be conveniently set by the user.


Efficiency

- Automation of time-consuming, labour-intensive manual operations
- ◆ Very large sample throughput
- Daily capacity can be accurately predicted

Functionality

- Automatic filling of up to four independently operating viscometers
- Automatic emptying and cleaning with up to two cleaning liquids
- Facility for mixed operation with two sample types dissolved in different media
- Operations controlled completely by PC
- Maximum safety in handling dangerous substances





Software


All PVS system configurations are operated via an interface by a conventional PC. The powerful and user-friendly PC program in its standard version also performs all necessary calculations for determining kinematic, dynamic, relative, reduced and inherent viscosity as well as the K-value, completely auto-matically on the basis of the measurements. Further substance characteristics can be obtained with additional software modules which can be interconnected to the basic software.

Software on Windows basis


Basic software


- Windows software, running on all conventional PCs and operating systems
- Parallel operation on up to eight places
- Calculation and presentation of
 - flow times and their average
 - standard deviation
 - kinetic energy correction e.g. Hagenbach correction or ISO 1628/6
 - absolute kinematic viscosity
 - absolute dynamic viscosity
 - relative viscosity
 - reduced viscosity (viscosity number)
 - inherent viscosity (logarithmic viscosity number)
 - K-value after Fickentscher

The program provides extensive user support. All parameters are input on the PC by mouse click and keyboard and are transferred to the PVS system via the interface.



Measurements can be printed in the form of a test protocol and stored in a file in ASCII format. The protocol covers all measurements of the day in chronological order and stores them in a file identified by date. The data can be viewed at any time, ensuring uninterrupted documentation. Further processing with other programs, such as MS Excel, and networking are readily available.






Compatibility with the laboratory environment

Good laboratory practice requires extensive independence of the measuring conditions by the operating staff. This means that as far as possible no critical measurement parameters such as for example sample descriptions and concentrations can be entered unmonitored. If inputs are required, these must be correspondingly authorised. The software in the PVS is capable of reading

in finished lists with all parameters provided by the LIMS; this also applies to concentration entries that can be read in directly from the communication-capable balances and dosing systems. The sample number can be transmitted to the display of the balance at the same time. This is effected from the measuring computer independently via a separate software module.


Access authorisations and documentation simplify the measurements

All users must log in with their own account and password and have limited access to the system dependent on a level of authorisation they enjoy. This means that shift staff can only read in finished lists and can activate a series of measurements with-out need for making any entries themselves. The automatic documentation records all results and users chronologically in daily protocol files and log files order as well as alterations made to parameters as required for example in FDA Standard 21CFR Part 11.



Example equipment

- Autosampler VAS 1/4 (complete)
- ♣ Analytic balances
- ♣ Burette 765
- LIMS software module (as per user specification)
- **◆ Stirrer block MRH 15**



Features

- Extensive support from GLP and 21 CFR part 11
- Automatic integration of the sample preparation
- Supports conventional balances and dosing systems
- User-specific connection to LIMS available


Solvent viscosity of plastics

Plastics industry today demands a large number of quality controls. From the raw material, through intermediate products and up to final processing, the average chain length of polymers with its decisive importance for quality, and its changes with

mechanical and thermal processing have to be checked again and again. The determination of solvent viscosity is here one of the most widely proven and sensitive methods.

Optimum system solutions for standard and special applications

The standard version of the PVS software already calculates automatically the relative and reduced viscosities, inherent viscosity and the K-value after Fickentscher. With the software module INV-DLL it is possible to determine the intrinsic viscosity by approximation formulae, either from one sample concentration (single-point method) or from different sample concentrations by extrapolation through linear regression. Intrinsic viscosity leads to mean molecular mass and from there to the chain length of a polymer. There are also taylor-made arrangements for polyolefines, such as polyethylene and polypropylene which can only be dissolved and tested at very high temperatures.



Single-point methods

Although these are based on approximate formulae with limited application and accuracy, they can be performed very rapidly on any PVS configuration since measurement at a single concentration only is required. The following methods are offered:

- · Point/slope method
- Schulz-Blaschke method
- · Huggins method
- · Solomon-Ciuta method
- Billmeyer method
- · Martin method
- · Maron method



Minimum specification


- **♣** Control unit PVS 1/1
- ♣ Measuring stand S 5
- ♣ KPG Ubbelohde viscometer
- ♣ PC with software modul INV-DLL
- ♣ Thermostat E 215 T with cover plate MD 15 V/K
- Through-flow cooler DLK 10 or cooling water connection

Almost all standardised methods based on solvent viscometry for the analysis of polymers can be realised in accordance with the specific application by means of the PVS. The most important standards for general solvent viscosity are ISO 1628/1, DIN 53728/1, ASTM D 2857; for PVC (polyvinyl chloride) are ISO 1628/2, DIN 53726, ASTM D 1243; for PE/PP (polyolefines)

are ISO 1628/3, DIN 53728/4, ASTM D 1601; for PC (polycarbonates) are ISO 1628/4, DIN 7744; for PET/PBT (polyester) are ISO 1628/5, DIN 53728/2, ASTM D 4603; for PMMA (polymethyl metacrylates) are ISO 1628/6; and for PA (polyamides) are ISO 307, DIN 53727, ASTM D 789.

Parallel linear regression

This is the fastest method on the market for performing linear regression. One measurement including cleaning takes only 25 minutes. The different sample concentrations required are evaluated virtually simultaneously on three to six places using independently operating measuring stands.




Minimum specification

- ♣ Control unit PVS 1/4
- ♣ 3 Measuring stands S 5
- ♣ 3 KPG Ubbelohde viscometers
- PC with software modul INV-DLL
- Thermostat PV 24 with cover plate D 20 V
- Through-flow cooler DLK 10 or cooling water connection

Serial linear regression

This method can be performed using only one place. Between the individual measurements the sample is successively diluted in the viscometer itself, using an automatic burette. A magnetic stirrer ensures rapid thermostating and uniformity of the sample after each concentration change. Apart from simplified operation, the main feature is highly precise dosing resulting in concentration series with excellent reproducibility.




Minimum specification

- ◆ Control unit PVS 1/1
- Measuring stand S 5
- Burette module BE and burette
- Ubbelohde dilution viscometer
- Magnetic stirrer
- PC with software module INV-DLL
- ♣ Thermostat E 215 T with cover plate MD 15 V/K
- Through-flow cooler DLK 10 or cooling water connection

Measuring system for polyolefines

This configuration permits simple and reliable measurement of the viscosity of polyethylene and polypropylene sulfide solutions at temperatures up to 160 °C. The sample, granulate or powder is placed directly into the dilution viscometer with integrated filters and is dissolved there, with the solvent required for this and for the subsequent dilution steps added precisely from a burette. Manual handling of hot samples or solvents, as usually required, has become unnecessary.



Typical specification

- ♣ Control unit PVS 1/4
- ♣ 4 Measuring stands S 5
- ♣ Burette module BE, 2 burettes
- 4 Ubbelohde dilution viscometers with integrated filters
- 2 Cleaning modules VRM 3D
- ◆ 2 temperable valve units UD 651
- ♣ 1 pump VRP
- PC with software module INV-DLL
- ♣ Thermostat PV 24 with cover plate D 20 V
- 4 Magnetic stirrers fitted in PV 24

Testing technical lubricants


The precision of the absolute kinematic viscosity measured using the LAUDA PVS system considerably exceeds the requirements determined in the standards. The LAUDA PVS in its optimised configured form is thus an ideal choice for the determination of viscosimetric characteristics such as of mi-

neral oils, other oils and their derivatives. The LAUDA PVS-System complies with many international standards, like kinematic viscometry DIN 51562/1-3, ASTM D 445-446, IP 71 and ISO 3104-3105 and for determination of the viscosity index ASTM D 2270 and ISO 2909.

Wide temperature ranges

Viscosity index

- ♣ Control unit PVS 1/8
- 8 Measuring stands S 5
- 8 Viscometers with aspirating tube (e.g. Cannon-Fenske routine)
- ♣ 4 Cleaning modules VRM 3
- ♣ PC with software module VID-DLL
- 2 Thermostats PV 24, for 40 and 100 °C each, plus cover plate D 20 V



Viscosity index

Mineral oil production involves regularly repeated measurements to determine and monitor viscosity and its variation at different temperatures. The software module VID-DLL provides a very efficient means for calculating the viscosity index. With one thermostat, the necessary measurements at 40 and 100 °C are performed after a temperature change. The method becomes particularly fast by using measuring stands in two thermostats where measurements at 40 and 100 °C are made almost simultaneously. With eight

places and four thermostats, a single procedure permits extremely convenient measurement of viscosities over a large temperature range of, for example, -20, 20, 40 and 100 °C.

Example of system for low-temperature viscosity down to -40 °C

- ♣ Control unit PVS 1/2
- 2 Measuring stands S 5
- ♣ Cleaning module VRM 3
- 2 Viscometers with aspirating tube (e.g. KPG Ubbelohde)
- ◆ PC
- ♣ Thermostat PVL 15 with cover plate D 15 V
- ♣ Through-flow cooler DLK 45 LiBus plus cold trap
- ◆ TEMP-DLL



Low-temperature viscosity

With the PVS system it is possible to determine the viscosity of fuel, e.g. for aircraft, under actual conditions of use to 60 °C. For the first time this is now also possible including automatic cleaning down to -20 °C. In addition to high-power LAUDA refrigeration units and specially insulated clear-view thermostats the me-


thod involves a cold trap for air drying. With software module TEMP-DLL the temperature of thermostats can be controlled and changed e.g. to determine the viscosity versus temperature dependency.

Integrated solvent recycling

In many cases only rinsing agents that are capable of completely dissolving polymers and oils can be used for the necessary thorough cleaning of viscometers. The rinsing solvents suitable for such a task are frequently chlorinated or toxic. Highly volatile rinsing solvents such as dichloromethane and

chloroform, as well as acetone, THF and some others can only be recovered automatically online by recooling, and this in the circuit directly out of waste bottle and back into the rinsing supply bottle.

Easy solvent handling



For this a PVS system with rinsing modules VRM 3D/HV is required in which the sample and the rinsing agent are aspirated by means of a vacuum membrane pump into a waste bottle. If this pump is fitted with a heat exchanger, the times between the rinsing processes can be utilised for recondensation of the solvent as the pump continues to run, reducing the pressure in the waste bottle down to the boiling pressure.

Evaporated solvent passes at the recooler where it is condensed and returned as purified liquid directly into the rinsing supply bottle: As the pump is used twice to aspirate the viscometer empty and to evaporate the rinsing agent, merely an additional recooler, and thermostating devices are required to cool the condensers and to keep the waste solvent boiling.

Apart from the money savings, handling of the solvent is also reduced as the necessity to refill is greatly lowered or is not even necessary, in case the sample solvent is the same as the rinsing solvent. Only LAUDA offers such an integrated solution that can be used for different equipment levels.

Example of a 4-place system including solvent recycling

- ♣ Control unit PVS 1/4
- ♣ 4 Measuring stands S 5
- 4 Ubbelohde with viscometer aspirating tube
- ♣ Thermostat PV 24 with coverplate D 20 V
- 2 Rinsing modules VRM 3D/HV
- Controlled vacuum pump with condensers
- ◆ PC
- ♣ Heating thermostat E 103
- ♣ Circulation chiller WKL 230

		Dichlormethane	Acetone	Chloroform	Ethylalcohol	
Boiling point	C°	40.0	56.5	61.0	78.5	
Temperature of cooler	C°	-10.0	-10.0	-10.0	-10.0	
Temperature in evaporation bath	C°	51.6	75.9	82.0	95.5	
Vapor pressure	hPa	350.0	350.0	230.0	200.0	
Evaporation time	min/100 ml	4.8	3.3	3.0	3.0	
Mass loss per cleaning	%	5.0	5.0	6.0	5.0	

Determining enzyme activity

Solutions of certain biological macromolecules alter their viscosity under the influence of enzymes as these effectively cut the dissolved molecular chains. Such situations permit very accurate evaluation of enzyme activity by measuring the variation of relative viscosity during the course of the enzyme reaction.

Automatic control of the measuring sequence

Apart from controlling the measurement sequence, the software module ENZ-DLL automatically calculates and outputs enzyme activity for hyaluronidase and cellulase from a comparison with reference measurements, in accordance with international pharmaceutical standards.

- Example of a system


 ◆ Control unit PVS 1/4
- ♣ 4 Measuring stands S 5
- ♣ 4 Ubbelohde viscometers
- with aspirating tube
- 2 Cleaning modules VRM 3PC with software module
- ENZ-DLL

 Thermostat PV 24 with
- cover plate D 20 V
- Through-flow cooler DLK 10 or cooling water connection

In addition the logarithmic relative viscosity against reaction time is presented graphically by the software module ENZ-DLL drawing the regression line through the test points and determin the reacting half-life.


Features

- Up to 99 individual measurements, cover of changes over a wide range of time
- Depending of system up to 8 measurements in parallel
- Reaction start (determined by mixing of solution) can be transferred by key stroke
- Setting of delay time and automatically start of measurements



Modules

PVS control units PVS 1		PVS 1/11/8
Max. number of places		18
PC interface		RS 232 C
Dimensions (WxDxH)	mm	340x270x105
Weight (net)	kg	4.6
Ambient temperature	°C	1045
Total loading	kW	0.1
Supply	V; Hz	90–240; 50/60



Measuring stands		S 5
Meniscus detection		optical (infrared)
Light detector control		digital (μΡ)
Sample temperature range	°C	65180*
Timing range	S	09999,99
Recommended flow timing		
range	S	301000
Viscosity range	mm²/s	0.350000
Timing resolution	S	0.01
Timing accuracy	ppm	1
Dimensions (WxDxH)	mm	90x90x500
Weight (net)	kg	4.5



^{*}higher temperatures to special order

Rinsing modules		VRM 3	VRM 3D/HV	VRM 3D/HT
Sample temperature range	°C	40100	40100	20165*
Viscosity range	mm²/s	0.3100**	0.31000**	0.3100**
Max. number of solvents		2	1	1
Dimensions (WxDxH)	mm	140x120x125	140x120x125	140x120x125
Weight (net)	kg	4.8	4.6	4.6








^{*}higher temperatures to special order

**can be ex	ktended using	samples	predilution	and	connection	set 8	3
-------------	---------------	---------	-------------	-----	------------	-------	---

Autosampler		VAS 1
Sample temperature range	°C	20135*
Viscosity range	mm²/s	0.3100
Syringe volume	°C	5
Max. number of samples (50 ml)		35/24
Max. number of samples (25 ml)		63/43
Max. number of places		2/4
Dimensions (WxDxH)	mm	1200x600x1200**
Weight (net)	kg	64**



^{*}thermostated syringe and sample stand with thermostat option
**complete system VAS 1/2 or VAS 1/4 with thermostat PV 24

Control unit PVS 1		Autosampler VAS 1 and accessories	Cat. No.
Туре	Cat. No.	Essential accessories	
PVS 1/1 incl. 1x ME 1 (for 1 place)	LMV 816	Coupling cap (1 x for each EG 066 necessary)	EZV 104
PVS 1/2 incl. 1x ME 2 (for 2 places)	LMV 810	Sealing rings (50 pieces) for EZV 104	EDF 124
PVS 1/4 incl. 2 x ME 2 (for 4 places)	LMV 813	Aluminium-plates (à 1000 pieces) for EG 066	EDF 092
PVS 1/6 incl. 3 x ME 2 (for 6 places)	LMV 813	Operating unit for VAS	EBE 038
PVS 1/8 incl. 4 x ME 2 (for 8 places)	LMV 815	Dosing syringe for VAS 1 (only spare part)	UD 442
1 VO 170 IIICI. 4 X W.L 2 (101 0 piaces)	LIN V 010	Dosing syringe to be used with LMVZ 958	UD 556
Each version includes		Heating block for dosing syringe	EBE 037
Windows software	LDVM 4014	Vacuum pump (controlled)	LMVZ 157
RS 232-cable for PC	EKS 037		
NO 202 00010 101 1 0	LIKO UUT	Software and accessories for sample preparation	
Plug-in extension cards		Туре	Cat. No.
2-place measurement module (ME 2) (for up to 2 places)	LMVZ 930	Analytic balance	EBK 006
Burette module (BE) (for up to 2 burettes 665)	LMVZ 932	Burette 765 (fully automatic)	EBK 003
Magnetic stirrer control module (for max. 4 stirrers) MRE		Stirring block MRH 15 (for 100 ml flasks)	EBK 010
Port control module (VAS/E)	LMVZ 943	Adapter rings for 50 ml flasks	EG 062
(in connection to the VAS autosampler)		Connecting cable for Burette 765	UK 253
(Software with own data base	LDVM 4022
		Software for customers data base, e.g. LIMS-systems.	LDVM 4023
Autosampler VAS 1 and accessories	Cat. No.		
Two-placesystem with Autosampler VAS 1/2	LMV 818	Measuring stands	Cat. No.
with software incl. PVS 1, 2 measuring stands S 5		Measuring stands S 5 (incl. cable and tubing)	LMVZ 948
cover plate for PV 24		PC and accessories	On request
Four-placesystem with Autosampler VAS 1/2	LMV 819		
incl. PVS 1, 4 measuring stands S 5		Software modules	
cover plate for PV 24		Software module INV-DLL (intrinsic viscosity)	LDVM 4015
		Software module VID-DLL (viscosity index to ISO 2909)	LDVM 4016
Essential accessories		Software module ENZ-DLL (enzyme activity)	LDVM 4017
Syringe Wash station 4 for VAS 1	LMVZ 970	Software module TEMP-DLL	LDVM 4023
(with one rinsing bottle)		(temperature control and dependence)	
Syringe Wash station 5 for VAS 1	LMVZ 971		
(with two rinsing bottles)		Dosing system	
Syringe Wash station 3 for VAS 1 (for hot solutions	LMVZ 972	Burette 765	EBK 003
with one rinsing bottles)		Burette set 1 for operation without VRM	LMVZ 931
Connecting set 3 for viscometers with aspiration tube	LMRZ 911	(for 1 burette and 1 viscometer)	
Connecting set 5 for viscometers without aspiration tube	LMRZ 914	Burette set 2 for operation with VRM	LMVZ 937
Connecting set 7 for automatic dilution	LMRZ 921	(for 1 burette and up to 2 viscometers)	
Sample rack PG 50, not heated	LMVZ 939		
(for 35 x 50 ml flasks, EG 062)	1 141/7 047	Additional accessories	
Sample rack PG 30, not heated	LMVZ 947	Draining rack (for filling and drying)	UU 004
(for 63 x 30 ml flasks)	LMV7.000	Bottle	LMVZ 934
Mounting set for stirrer block MRH 15	LMVZ 969	Connecting cap, small, silicone	HKA 001
Sample rack PGH 24 (heated, max.160°C)	LMVZ 946	Connecting cap, large, silicone	HKA 002
(for 24 x 50 ml flasks, EG 062)	EG 062	Connecting cap, small, viton	HKA 147
Flasks (50 ml), with GL32 thread for PG 50		Connecting cap, large, viton	HKA 148
Filter element for EG 062 (1 for each EG 062) Coupling cap (1 x for each EG 062 necessary)	LMVZ 958 EZV 100	Connecting cap, large, silicone (for dilution series)	HKA 118
Sealing rings (50 pieces) for EZV 100	EDF 122	Silicon tubes, 3 x 1.5 mm	RKJ 014
Aluminium-plates (à 1000 pieces) for EG 062	EDF 122	Viton tube, 3 x 1.5 mm (for sulfuric acid)	RICJ 020 UK 237
Flasks (30 ml), with GL32 thread for PG 30	EG 066	Connecting cable burette <-> PVS	
i lasks (50 iiii), with OLSZ tilleau ioi FO 50	_0 000	PTFE stirrer	EZ 195

Automatic cleaning


VRM 3	
Туре	Cat. No.
VRM 3 set (for up to 2 measuring stands	LMR 909
or syringe wash stations)	
Standard parts	
Cleaning module VRM 3	UD 640
Tubing set 2	LMRZ 903
2 connecting caps for glass bottles (GL 45)	LMRZ 907
Connecting cable VRM <-> PVS	UK 230
Further accessories	Cat. No.
Glass funnel for liquid samples	EG 060
Funnel for granulate	HX 488
Filter insert for filter EG 060	UD 410
Filter for solvent	UD 404
Glass bottle 1000 ml (GL 45)*	EG 058
Glass bottle 2000 ml (GL 45)**	EG 059
Glass bottle 5000 ml (GL 45)*	EG 6064
Connecting caps for glass bottles (GL 45)	LMRZ 907
One way filling syringe (100 pcs)	LMRZ 918
Filter for suction	UD 513
*recommended for solvent **recommended for waste	

Fixing sets for VRM modules*	Cat. No.
for PV 15 / PVL 15	LMRZ 904
for PV 24 / PVL 24	LMRZ 905
for PV 36	LMRZ 906
for E 215 T* (for 1 VRM)	LMRZ 908
for E 215 T* (for 2 VRM)	LMRZ 916

^{*}for mounting on a thermostat $% \left(1\right) =\left(1\right) \left(1\right) \left($

Type Cat. No. VRM 3D set (for max. two measuring stands or syringe wash stations) Standard parts Cleaning module VRM 3D UD 650 Valve unit HV UD 652 (for particle containing or high viscous samples) Valve unit HT (for hot samples) UD 651 Tubing set 2 LMRZ 903 Connecting caps for glass bottles (GL 45) LMRZ 907 Connecting cable VRM <-> PVS UK 230 Essential accessories External vacuum pump VRP (controlled) LMRZ 809 for up to 3 VRM 3D External vacuum pump (not controlled) LMVZ 135 for up to 3 VRM 3D Glass bottle 2000 ml (GL 45) EG 059 Connecting set 1 (only for viscometers LMRZ 909 with aspirating tube) Connecting set 6 for autom. dilution LMRZ 919 Connecting set 8 for high viscous samples LMRZ 922		
VRM 3D set (for max. two measuring stands or syringe wash stations) Standard parts Cleaning module VRM 3D UD 650 Valve unit HV UD 652 (for particle containing or high viscous samples) Valve unit HT (for hot samples) UD 651 Tubing set 2 LMRZ 903 Connecting caps for glass bottles (GL 45) LMRZ 907 Connecting cable VRM <-> PVS UK 230 Essential accessories External vacuum pump VRP (controlled) LMRZ 809 for up to 3 VRM 3D External vacuum pump (not controlled) LMVZ 135 for up to 3 VRM 3D Glass bottle 2000 ml (GL 45) EG 059 Connection sets Cat. No. Connecting set 1 (only for viscometers LMRZ 909 with aspirating tube) Connecting set 6 for autom. dilution LMRZ 919	VRM 3D/HT ·VRM 3D/HV	
Standard parts Cleaning module VRM 3D Valve unit HV UD 652 (for particle containing or high viscous samples) Valve unit HT (for hot samples) UD 651 Tubing set 2 Connecting caps for glass bottles (GL 45) Connecting cable VRM <-> PVS UK 230 Essential accessories External vacuum pump VRP (controlled) for up to 3 VRM 3D External vacuum pump (not controlled) for up to 3 VRM 3D Glass bottle 2000 ml (GL 45) Connection sets Cat. No. Connecting set 1 (only for viscometers with aspirating tube) Connecting set 6 for autom. dilution LMRZ 919	Туре	Cat. No.
Standard parts Cleaning module VRM 3D Valve unit HV (for particle containing or high viscous samples) Valve unit HT (for hot samples) UD 651 Tubing set 2 Connecting caps for glass bottles (GL 45) Connecting cable VRM <-> PVS UK 230 Essential accessories External vacuum pump VRP (controlled) for up to 3 VRM 3D External vacuum pump (not controlled) for up to 3 VRM 3D Glass bottle 2000 ml (GL 45) Connection sets Cat. No. Connecting set 1 (only for viscometers with aspirating tube) Connecting set 2 (only for viscometers without aspirating tube) Connecting set 6 for autom. dilution LMRZ 919	VRM 3D set (for max. two measuring stands	LMR 907
Cleaning module VRM 3D Valve unit HV (for particle containing or high viscous samples) Valve unit HT (for hot samples) UD 651 Tubing set 2 Connecting caps for glass bottles (GL 45) Connecting cable VRM <-> PVS UK 230 Essential accessories External vacuum pump VRP (controlled) for up to 3 VRM 3D External vacuum pump (not controlled) for up to 3 VRM 3D Glass bottle 2000 ml (GL 45) Connection sets Cat. No. Connecting set 1 (only for viscometers with aspirating tube) Connecting set 2 (only for viscometers without aspirating tube) Connecting set 6 for autom. dilution LMRZ 919	or syringe wash stations)	
Cleaning module VRM 3D Valve unit HV (for particle containing or high viscous samples) Valve unit HT (for hot samples) UD 651 Tubing set 2 Connecting caps for glass bottles (GL 45) Connecting cable VRM <-> PVS UK 230 Essential accessories External vacuum pump VRP (controlled) for up to 3 VRM 3D External vacuum pump (not controlled) for up to 3 VRM 3D Glass bottle 2000 ml (GL 45) Connection sets Cat. No. Connecting set 1 (only for viscometers with aspirating tube) Connecting set 2 (only for viscometers without aspirating tube) Connecting set 6 for autom. dilution LMRZ 919		
Valve unit HV (for particle containing or high viscous samples) Valve unit HT (for hot samples) UD 651 Tubing set 2 Connecting caps for glass bottles (GL 45) Connecting cable VRM <-> PVS UK 230 Essential accessories External vacuum pump VRP (controlled) for up to 3 VRM 3D External vacuum pump (not controlled) for up to 3 VRM 3D Glass bottle 2000 ml (GL 45) Connection sets Cat. No. Connecting set 1 (only for viscometers with aspirating tube) Connecting set 2 (only for viscometers without aspirating tube) Connecting set 6 for autom. dilution LMRZ 919		
(for particle containing or high viscous samples) Valve unit HT (for hot samples) Tubing set 2 Connecting caps for glass bottles (GL 45) Connecting cable VRM <-> PVS UK 230 Essential accessories External vacuum pump VRP (controlled) for up to 3 VRM 3D External vacuum pump (not controlled) for up to 3 VRM 3D Glass bottle 2000 ml (GL 45) Connection sets Cat. No. Connecting set 1 (only for viscometers with aspirating tube) Connecting set 2 (only for viscometers without aspirating tube) Connecting set 6 for autom. dilution LMRZ 919	Cleaning module VRM 3D	UD 650
Valve unit HT (for hot samples) Tubing set 2 Connecting caps for glass bottles (GL 45) Connecting cable VRM <-> PVS UK 230 Essential accessories External vacuum pump VRP (controlled) for up to 3 VRM 3D External vacuum pump (not controlled) for up to 3 VRM 3D Glass bottle 2000 ml (GL 45) Connection sets Cat. No. Connecting set 1 (only for viscometers with aspirating tube) Connecting set 2 (only for viscometers without aspirating tube) Connecting set 6 for autom. dilution LMRZ 919	Valve unit HV	UD 652
Tubing set 2 Connecting caps for glass bottles (GL 45) Connecting cable VRM <-> PVS UK 230 Essential accessories External vacuum pump VRP (controlled) for up to 3 VRM 3D External vacuum pump (not controlled) for up to 3 VRM 3D Glass bottle 2000 ml (GL 45) Connection sets Cat. No. Connecting set 1 (only for viscometers with aspirating tube) Connecting set 2 (only for viscometers without aspirating tube) Connecting set 6 for autom. dilution LMRZ 919	(for particle containing or high viscous samples)	
Connecting caps for glass bottles (GL 45) Connecting cable VRM <-> PVS UK 230 Essential accessories External vacuum pump VRP (controlled) for up to 3 VRM 3D External vacuum pump (not controlled) for up to 3 VRM 3D Glass bottle 2000 ml (GL 45) Connection sets Cat. No. Connecting set 1 (only for viscometers with aspirating tube) Connecting set 2 (only for viscometers without aspirating tube) Connecting set 6 for autom. dilution LMRZ 919	Valve unit HT (for hot samples)	UD 651
Connecting cable VRM <-> PVS Essential accessories External vacuum pump VRP (controlled) for up to 3 VRM 3D External vacuum pump (not controlled) for up to 3 VRM 3D Glass bottle 2000 ml (GL 45) Connection sets Cat. No. Connecting set 1 (only for viscometers with aspirating tube) Connecting set 2 (only for viscometers without aspirating tube) Connecting set 6 for autom. dilution LMRZ 919	Tubing set 2	LMRZ 903
Essential accessories External vacuum pump VRP (controlled) for up to 3 VRM 3D External vacuum pump (not controlled) for up to 3 VRM 3D Glass bottle 2000 ml (GL 45) Connection sets Connecting set 1 (only for viscometers with aspirating tube) Connecting set 2 (only for viscometers without aspirating tube) Connecting set 6 for autom. dilution LMRZ 919	Connecting caps for glass bottles (GL 45)	LMRZ 907
External vacuum pump VRP (controlled) for up to 3 VRM 3D External vacuum pump (not controlled) for up to 3 VRM 3D Glass bottle 2000 ml (GL 45) Connection sets Cat. No. Connecting set 1 (only for viscometers with aspirating tube) Connecting set 2 (only for viscometers without aspirating tube) Connecting set 6 for autom. dilution LMRZ 919	Connecting cable VRM <-> PVS	UK 230
External vacuum pump VRP (controlled) for up to 3 VRM 3D External vacuum pump (not controlled) for up to 3 VRM 3D Glass bottle 2000 ml (GL 45) Connection sets Cat. No. Connecting set 1 (only for viscometers with aspirating tube) Connecting set 2 (only for viscometers without aspirating tube) Connecting set 6 for autom. dilution LMRZ 919		
for up to 3 VRM 3D External vacuum pump (not controlled) for up to 3 VRM 3D Glass bottle 2000 ml (GL 45) Connection sets Cat. No. Connecting set 1 (only for viscometers with aspirating tube) Connecting set 2 (only for viscometers without aspirating tube) Connecting set 6 for autom. dilution LMRZ 919	Essential accessories	
External vacuum pump (not controlled) for up to 3 VRM 3D Glass bottle 2000 ml (GL 45) Connection sets Cat. No. Connecting set 1 (only for viscometers LMRZ 909 with aspirating tube) Connecting set 2 (only for viscometers LMRZ 910 without aspirating tube) Connecting set 6 for autom. dilution LMRZ 919	External vacuum pump VRP (controlled)	LMRZ 809
for up to 3 VRM 3D Glass bottle 2000 ml (GL 45) Connection sets Connecting set 1 (only for viscometers with aspirating tube) Connecting set 2 (only for viscometers without aspirating tube) Connecting set 6 for autom. dilution LMRZ 919	for up to 3 VRM 3D	
Connection sets Connecting set 1 (only for viscometers with aspirating tube) Connecting set 2 (only for viscometers without aspirating tube) Connecting set 6 for autom. dilution EG 059 Cat. No. LMRZ 909 WHRZ 910 LMRZ 910 LMRZ 910	External vacuum pump (not controlled)	LMVZ 135
Connection sets Connecting set 1 (only for viscometers with aspirating tube) Connecting set 2 (only for viscometers without aspirating tube) Connecting set 6 for autom. dilution Cat. No. LMRZ 909 LMRZ 919	for up to 3 VRM 3D	
Connecting set 1 (only for viscometers with aspirating tube) Connecting set 2 (only for viscometers without aspirating tube) Connecting set 6 for autom. dilution LMRZ 919	Glass bottle 2000 ml (GL 45)	EG 059
Connecting set 1 (only for viscometers with aspirating tube) Connecting set 2 (only for viscometers without aspirating tube) Connecting set 6 for autom. dilution LMRZ 919		
Connecting set 1 (only for viscometers with aspirating tube) Connecting set 2 (only for viscometers without aspirating tube) Connecting set 6 for autom. dilution LMRZ 919		
with aspirating tube) Connecting set 2 (only for viscometers without aspirating tube) Connecting set 6 for autom. dilution LMRZ 919	Connection sets	Cat. No.
Connecting set 2 (only for viscometers LMRZ 910 without aspirating tube) Connecting set 6 for autom. dilution LMRZ 919		LMRZ 909
without aspirating tube) Connecting set 6 for autom. dilution LMRZ 919	with aspirating tube)	
Connecting set 6 for autom. dilution LMRZ 919	Connecting set 2 (only for viscometers	LMRZ 910
-		
Connecting set 8 for high viscous samples LMRZ 922	Connecting set 6 for autom. dilution	LMRZ 919
	Connecting set 8 for high viscous samples	LMRZ 922

Viscometers for offline cleaning



Ubbelohde viscometers

ISO 3105, DIN 51562, BS 188, NFT 60-100.

Filling volume: 15...20 ml Total length: 290 mm approx.

Accuracy: \pm 0.1%, calibrated for absolute measu-

rement, for automatic measurement. Also available in ASTM version.

Туре	K	Cat. No.
0	0.001	EGV 690
0c	0.003	EGV 700
0a	0.005	EGV 701
1	0.01	EGV 702
Ic	0.03	EGV 703
la	0.05	EGV 691
II	0.1	EGV 704
IIc	0.3	EGV 705
lla	0.5	EGV 692
III	1	EGV 706
IIIc	3	EGV 707
Illa	5	EGV 693
IV	10	EGV 708
IVc	30	EGV 699
IV a	50	EGV 698




Cannon-Fenske-Routine viscometers

ISO 3105, ASTM D 2515, BS 188

Filling volume: 5...10 ml Total length: 245 mm approx.

Accuracy: ± 0.2%, calibrated for absolute measurement, for automatic measurement.

Туре	K	Cat. No.
25	0.002	EGV 860
75	0.008	EGV 861
50	0.004	EGV 862
100	0.015	EGV 863
150	0.035	EGV 864
250	0.1	EGV 865
350	0.5	EGV 866
300	0.25	EGV 867
400	1.2	EGV 878
450	2.5	EGV 869
500	8	EGV 870
600	20	EGV 871




Micro-Ubbelohde viscometers

Filling volume: 2...3 ml Total length: 290 mm approx.

Accuracy: ± 0.5%, calibrated for absolute measurement, for automatic measurement.

Туре	K	Cat. No.
1	0.01	EGV 718
Ic	0.03	EGV 719
II	0.1	EGV 720
IIc	0.3	EGV 721
III	1	EGV 722



Micro-Ostwald viscometers

Recommended with pronounced foaming

and/or small liquid quantities

Filling volume: 2 ml

Total length: 290 mm approx.


Accuracy: ± 0.2%, calibrated for absolute measurement, for automatic measurement.

Type	K	Cat. No.
T	0.002	EGV 820
Ic	0.008	EGV 821
II	0.004	EGV 822
IIc	0.015	EGV 823
III	0.035	EGV 824

Viscometer holders and accessories:	
Туре	Cat. No.
Adapter MUO, required for use of Micro-Ubbelohde or Micro-Ostwald viscometers	HBK 532
Ubbelohde viscometer holder*	UG 003
Cannon-Fenske viscometer holder*	UG 084
Micro-Ostwald viscometer holder*	UG 094

^{*}holder fits only LAUDA viscometers

Viscometers for offline cleaning



Dilution-Ubbelohde viscometers

For determining the intrinsic viscosity

Filling volume: 15...75 ml Total length: 290 mm approx. Accuracy: ± 0.1%, uncalibrated, for

automatic measurement.

Type	K	Cat. No.
0	0.0011	EGV 920
0a	0.005	EGV 921
0c	0.001	EGV 922
1	0.01	EGV 923*
lc	0.01	EGV 924
II	0.1	EGV 925


*also available with integrated filter

Calculation formula:

Kinematic viscosity =

viscometer constant (K) x flow time*
*not taken into account Hagenbach correction

Viscometers for online cleaning



Ubbelohde viscometers

ISO 3105, DIN 51562, BS 188,

NFT 60-100.

With filling and cleaning tube Filling volume: 18...22 ml Total length: 290 mm approx.

Accuracy: ± 0.1%, calibrated for absolute measu-

rement, for automatic measurement.
Only in conjunction with VRM modules.

Also available in ASTM version.

71.00 a.caa		
Туре	K	Cat. No.
25	0.002	EGV 930
75	0.008	EGV 931
50	0.004	EGV 932
100	0.015	EGV 933
150	0.035	EGV 934
250	0.1	EGV 935
350	0.5	EGV 936
300	0.25	EGV 937
400	1.2	EGV 938

Cannon-Fenske-Routine viscometers


ISO 3105, DIN 51562, BS 188. With filling and $\,$

cleaning tube

Filling volume: 5..10 ml Total length: 245 mm approx.

Accuracy: ± 0.2%, calibrated for absolute measu-

rement, for automatic measurement.
Only in conjunction with VRM modules.



Type	K	Cat. No.
25	0.002	EGV 950
75	0.008	EGV 951
50	0.004	EGV 952
100	0.015	EGV 953
150	0.035	EGV 954
250	0.1	EGV 955
350	0.5	EGV 956
300	0.25	EGV 957
400	1.2	EGV 958
450	2.5	EGV 959
500	8	EGV 960
600	20	EGV 861

Clear-view thermostats

LAUDA clear-view thermostats are available in various bath sizes and for various bath depths, as both a Master and a Command version. All clear-view thermostats are equipped with a Varioflex pump. The larger cooling coils, which come

built in as standard, offer a more effective cooling capacity, especially for large baths. An easily-accessible drain tap on the side of the thermostat plus handle simplify mobility.

LAUDA Proline clear-view thermostats ensure accurate and reliable thermostating of objects placed in the bath. Because of the transparent front panel or the transparent bath vessel they are particularly suitable for viscometry. Provision of a special cover pla-

te makes the thermostats suitable for use with capillary viscometers and the stands of the PVS viscosity measuring system. In addition the models series PV and PVL incorporate an unique 2-chamber system; separation into a measurement and a thermostating chamber offers decisive advantages in the

measurement chamber: constant liquid level, very small temperature gradient, maximum temperature stability. In addition they incorporate a powerful pump and connectors for connecting to LAUDA through-flow coolers. The Ecoline Staredition thermostats E 115 T and E 215 T offer a lower-priced alternative in the temperature range up to 100 °C.



		100 °C	with RS 232	230 °C	24 litre	36 litre	-60 °C	-60 °C
Thermostat		E 115 T	E 215 T	PV 15	PV 24	PV 36	PVL 15	PVL 24
Working temp. range	°C	20100	20100	30230	30230	30230	30100	30100
Operating temp. range	°C	-20100	-20100	0230	0230	0230	-60100	-60100
Temperature stability	± °C	0.02	0.01	0.01	0.01	0.01	0.01	0.01
Resolution of setting	°C	0.1	0.1/0.01	0.1/0.01	0.1/0.01	0.1/0.01	0.1/0.01	0.1/0.01
Resolution of indication	°C	0.01	0.05	0.01	0.01	0.01	0.01	0.01
Safety fittings*		III, FL	III, FL	III, FL	III, FL	III, FL	III, FL	III, FL
Heater power	kW	1.5	2.25	3.5	3.5	3.5	3.5	3.5
Pump								
Туре	_				 pressure pur 	•		
Pump pressure max.	bar	0.4	0.4	0.8	0.8	0.8	0.8	0.8
Pump flow max. (pressure)	l/min	17	17	25	25	25	25	25
Bath								
Volume	1	1015	1015	1115	1924	2836	1115	1924
Opening (WxD)	mm	275x130	275x130	230x135	405x135	585x135	230x135	405x135
Bath depth	mm	310	310	320	320	320	320	320
Usable depth	mm	290	290	285	285	285	285	285
Material	mm	I polycart	oonateI		stair	nless steel/glass		I
Base area (WxD)	mm	428x142	428x142	506x282	740x282	1040x282	506x282	740x282
Cat. No.		LCD 0263	LCD 0264	LCD 0276	LCD 0278	LCD0280	LCD 0282	LCD 0284

^{*}III, FL for use with flammable and non-flammable liquids

For all clear-view thermostats there are various cover plates whose use is strongly recommended.

				Use	with	Together with
Туре	Cat. No.	No. meas. places	No. therm. places*	PVS without VRM	PVS with VRM	Comb. No.
MD 15 V	LCZ 041	2	0	•	•	1, 2, 3
MD 15 V/K	LCZ 040	1	2	•		1, 2, 3
D 15 V	LTZ 045	2	0	•	•	4, 5, 6, 13
D 15 V/K	LTZ 048	1	2	•		4, 5, 6, 13
D 20 V	LTZ 046	4	0	•	•	7, 8, 9, 14, 16
D 20 V/K	LTZ 049	3	3	•		7, 8, 9, 14, 16
D 30 V	LTZ 047	6	0	•	•	10, 11, 12

^{*}Thermostating places can be used for preliminary thermostating of filled viscometers in order to save time during measurement. They are only useful if the viscometer is changed for each measurement.

Through-flow coolers

LAUDA through-flow coolers serve as additional coolers in order to cool down heating thermostats or any baths to temperatures below ambient temperature. Through-flow coolers are used mainly for thermostats and are integrated into the cooling circuit. The high-power model DLK 45 LiBus even reaches low temperatures down to -40 $^{\circ}$ C.

Through-flow coolers not only render mains water cooling unnecessary, they also prevent undesirable fluctuations in the flow rate and ensure constant cooling water temperature. This leads to optimum temperature control over the entire measurement period which in turn has a very positive effect on the accuracy and reproducibility of the measurements.



		-15 °C	-30 °C	-40 °C
Through flow coolers		DLK 10	DLK 25	DLK 45 LiBus
Working temperature range	°C	-15150	-30150	-40150
Ambient temperature range	°C	540	540	540
Cooling output at, gross 20 °C	kW	0.25	0.33	1.10
0 °C	kW	0.20	0.28	0.95
-10 °C	kW	0.10	0.25	0.85
-20 °C	kW		0.22	0.75
-30 °C	kW		0.20	0.55
-40 °C	kW			0.30
Heat exchanger connection		M	16 x 1.13 mm dia. nipple	
Special functions		control	connection for 230 V 50/60 Hz	
Overall dimensions (WxDxH)	mm	200x400x320	290x540x330	470x560x430
Weight	kg	17	33	63
Supply*	V; Hz	230; 50/60	230; 50	230; 50
Power loading	kW	0.2	0.5	0.9
CatNr.		LFD 010	LFD 108	LFD 111

^{*}Protection Class 1 to VDE 0106

Magnetic stirrer sets for determining the intrinsic viscosity by serial regression

			No. dilution	No. magnetic	
Туре	Thermostat	Meas. places used	viscometers	stirrer sets	Cat. No.
1-place set	E 115 T, E 215 T	1-2	1-2	1-2	LMZ 841
2-place fitting set*	PV 15	1-2	1-2	1	LMVZ 967
4-place fitting set*	PV 24	1-4	1-4	1	LMVZ 968

^{*}built into thermostat, cannot be retrofitted

Туре	Cat. No.
Cold Trap	LMRZ 915
(necessary for air-drying at T < 0 °C	
for connecting to a PVS-System)	
Tubes for connecting thermostat/DLK	
Silicon tube 8 mm i. dia. (9 mm insulated)	LZS 001
Silicon tube 11 mm i. dia. (9 mm insulated)	LZS 007

Essential accessories:

Туре	Cat. No.	Comb.	
Connecting cabel	UK 263	3, 6, 9, 12 16	

Recommended accessories:

Туре	Cat. No.	Comb.
Silicone tubing* (per m,		
min. recommended: 2 m)	LZS 007	3. 6. 9. 12 16

^{*}not when using Silicone oil



The table below shows the system combinations for thermostating the viscometers:

				Clear-view thermostats/external baths		Cooling	
Comb. No.	Tmax (°C)	Tmin (°C)	Max. No. places	Type	Cat. No.	Туре	Cat. No.
1	100	30	2	E 115 T	LCD 0263		
				E 215 T	LCD 0264		
2	100	20	2	E 115 T	LCD 0263	Tap water	
				E 215 T	LCD 0264	Tap water	
3	100	5	2	E 115 T	LCD 0263	DLK 10	LFD 010
				E 215 T	LCD 0264	DLK 10	LFD 010
4	230	30	2	PV 15	LCD 0276		
5	230	25	2	PV 15	LCD 0276	Tap water	
6	230	10	2	PV 15	LCD 0276	DLK 10	LFD 010
7	230	30	4	PV 24	LCD 0278		
8	230	25	4	PV 24	LCD 0278	Tap water	
9	230	15	4	PV 24	LCD 0278	DLK 10	LFD 010
10	230	30	6	PV 36	LCD 0280		
11	230	25	6	PV 36	LCD 0280	Tap water	
12	230	15	6	PV 36	LCD 0280	DLK 25	LFD 108
13	100	-20	2	PVL 15	LCD 0282	DLK 25	LFD 108
14	100	-20	4	PVL 24	LCD 0284	DLK 25	LFD 108
15	100	-40	2	PVL 15	LCD 0282	DLK 45 LiBus	LFD 111
16	100	-40	4	PVL 24	LCD 0284	DLK 45 LiBus	LFD 111

- 1. All values Tmax and Tmin are based on an ambient temperature of 20 $^{\circ}\text{C}.$
- 2. All values for tap water cooling are based on a water temperature of maximal 18 $^{\circ}\text{C}.$
- 3. The temperature range Tmin...Tmax usually requires two different bath liquids.

Background illumination and accessories

Туре	CatNr.
AL Atherman lamp 15 (for PV 15 and PVL 15)	LTZ 001
AL Atherman lamp 20 (for PV 24 and PVL 24)	LTZ 002
AL Atherman lamp 30 (for PV 36)	LTZ 003
Filter for thermostating bath	EG 065
Flow indicator (necessary when using EG 065)	EZ 204
Viscometer holder for 2-legged capillaries	EZ 054
(for manual measurement only)	

21 CFR-11

The guideline 21 CRF, part 11, issued by the American FDA authority, regulates the technical and organisational requirements which must be fulfilled in order to use electronic data and documents instead of paper in the development, authorisation and production processes.

Billmeyer formula (IV value according to Billmeyer)

Serves the purpose of the approximate calculation of the intrinsic viscosity of polyesters and others. No additional polymer characteristic data required.

$$v_{\text{int}} = \frac{1}{4} v_{\text{red}} + \frac{3x \ln v_{\text{rel}}}{4C}$$

Dynamic viscosity

Is essential to the shear flows, viscosity coefficient, between shearing strength τ and speed gradient D in τ = η D, and has the unit mPas (formerly centipoise, cps).

FDA

Abbreviation for the Food and Drug Administration, United States. Prescribes binding guidelines for the development and production of phar-maceutical products: is internationally valid.

Glass viscometer

Viscometers made from glass, bearing different various designs, stan-dardised in ISO 3105. The most common for automatic measurements is the Ubbelohde version with ventilation pipe.

GLP

Abbreviation for "Good Laboratory Practice". Specifications initiated by the FDA American authority for laboratories and producers (e.g. of pharmaceuticals) regarding how tests and measurements are to be cleanly planned, performed and monitored. The guidelines have a legal character in many countries.

Hagen Poiseuill's Law

(fundamental equation of the capillary viscometry)

This forms the basis for viscometry in capillary viscometers. If the differential pressure is generated by a height difference in front of and behind the capillaries, hence: $\mathbf{v}_{\text{kin}} = \mathbf{k} \ \mathbf{x} \ t$, (k: capillary constant, t: measured run of a defined volume of liquid). In the case of very short times, the non-dissipated kinetic energy must be taken into consideration (kinetic energy / Hagenbach correction)

Huggins formula

(IV value according to Huggins)

Serves the purpose of the approximate calculation of the intrinsic viscosity of polystyrenes and others, for example. $K_{\rm H}$ is an additional constant dependent on polymer.

$$v_{\text{int}} = \frac{\sqrt{1 + 4K_{\text{H}} \times v_{\text{spec}}} - 1}{2xCxK_{\text{H}}}$$

Inherent viscosity

(logarithmic viscosity coefficient)

Is the natural logarithm of the relative viscosity based on concentration C of the dissolved substance $v_{\rm red} = :\ln v_{\rm rel} / C;$ unit: cm³/g = 100 dl/g

Intrinsic viscosity

(limiting viscosity number, Staudinger index, IV value)

Is the limiting value of the reduced/inherent viscosity for the case of infinitely severely diluted solvents at disappearing shearing strengths:

$$v_{\text{int}} = \lim_{C \to 0} v_{\text{red}} / v_{\text{int}} = \lim_{C \to 0} v_{\text{inh}}$$

It is determined by measuring the $v_{\rm nd}$ as a function of the concentration and extrapolation on C=0. For many polymers, there are approximation conditions based on the measurement of only one concentration usually specified in standards.

Kinematic viscometry

Describes the quotients of the dynamic viscosity by the density: $\nu_{\mbox{\tiny kin}} = \eta / \rho$ and has the unit mm²/s (formerly: centistokes, cst).

Kinetic energy correction: (Hagenbach correction)

If, in the case of short rundown times, there is a necessary correction of the Hagenbach-Poiseuill's Law, and it takes into consideration the kinetic energy not converted into friction warmth in a capillary viscometer.

Corrected viscosity
$$v = k x (t - \Delta t)$$

Correction factor
$$\Delta t = \frac{E}{k \times t^2}$$

Correction factor ISO 1628/6:
$$\Delta t = k/t - v_{ref}(t)$$

K value (according to Fickentscher)

A traditionally-used relative mass for the mole masse for PVC and PVA.

$$K = \frac{a-1+\sqrt{1+\left(\frac{2}{C}+2+a\right)xa}}{150+300C}$$

with: $a = 1.5 \times \log v_{rol}$

LIMS

Abbreviation for Laboratory Management System. Describes a system for the control and management of laboratory data, determined by various measuring devices.

Mark-Houwink formula

Provides the relation between medium mole masse (weight means) of the dissolved polymer chains and the intrinsic viscosity. For the absolute mole mass, the proportional constant K and the exponent a can be entered. These depend on the polymer and the solvent, and can be taken from the literature.

$$M = \left(\frac{v_{int}}{k}\right)^{\frac{1}{a}}$$

Martin formula (IV value according to Martin)

Serves the purpose of the approximate calculation of the intrinsic viscosity of celluloses and others, e. g. K is an additional constant dependent on polymer.

$$\log v_{red} = \log v_{int} + k \times v_{int} \times C$$

Reduced viscosity (viscosity coefficient)

Is the specific viscosity based on the concentration C of the dissolved substance.

$$v_{red} = : v_{sp} / C; unit: cm^3/g = 100 dl/g$$

Relative viscosity

Is the ratio of the dynamic viscosity η of the solvent to that of the solvent η_s . In the case of severely diluted solvents, this corresponds almost to the ratio of the kinematic viscosities:

$$V_{rel} = \eta / \eta_s \approx \nu / \nu_s$$

Schulz-Blaschke

(IV value according to Schulz-Blaschke)

Serves the purpose of the approximate calculation of the intrinsic viscosity of celluloses, polyolefines and others. K₁ is an additional constant dependent on polymer.

$$v_{\text{int}} = \frac{v_{\text{red}}}{1 + K_1 \times C \times v_{\text{red}}}$$

Solomon-Ciuta formula

(Solomon-Ciuta)

Serves the purpose of the approximate calculation of the intrinsic viscosity of PMMA and others. No additional polymer characteristic data required.

$$v_{\text{int}} = \frac{\sqrt{2 \, \text{x} \, (v_{\text{red}} \, \text{x} \, \text{C} - 1 \, \text{n} \, v_{\text{rel}})}}{C}$$

Specific viscosity (relative viscosity increase)

Is the relative viscosity minus one:

$$V_{sp} = V_{rel} - 1$$

Viscosity

Properties of a substance (in this case: of liquids) to flow and become irreversibly deformed under the influence of a stress. Flow energy is converted into warmth.

Viscosity index (for mineral oil products)

Is calculated from the viscosities measured at two different temperatures (40 and 100 °C) according to ISO 2909 and ASTM D 2270. Is a standard for the thermal behaviour of various oils. The higher the viscosity index of an oil is, the less it changes its viscosity at various temperatures.

Our product lines:

Thermostats \cdot Circulation chillers \cdot Water baths Heating and Cooling systems

Viscometers · Tensiometers

