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*Milli-Q® IQ 7003 및 Milli-Q® IQ 7005 시스템.

Compliance with Quality Assurance demands.

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기관	순수 품질/ 등급	초순수 품질 / 등급
유럽 약전	Purified water	Highly purified water
미국 약전	Purified water	Purified water in bulk
일본 약전	Purified water	Purified water
중국 약전	Purified water	Purified water
ASTM® D1193	Type II water	Type I water, Grade B
ISO® 3696	Grade 2 water	Grade 1 water
Chinese National Standard GB/T 6682	Level 2 water	Grade 1 water
Chinese National Standard GB/T 33087	N/A	Ultrapure water for instrumental analysis
JIS K 0557	A2 water	A4 water
Clinical and Laboratory Standards Institute® (CLSI®)	N/A	Clinical Laboratory Reagent Water (CLRW)

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
Technical Appendix



Milli-Q® IQ 7003/05/10/15 Water Purification Systems

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Production flow schematic


Distribution flow schematic


Purification unit and storage tank specifications

Purification unit

Storage tank

Tubing and port requirements

Parameter	Description
Feed water connection	½" Gaz
Distance from feed water port	Maximum 5 m
Distance from purification unit to POD dispenser	Maximum 5 m
Distance from purification unit to tank	Maximum 5 m
Power entry	Connection IEC 13
ON/OFF switch	Available on the unit
Water sensor port	Maximum 3.3 VDC
Tank level adapter port	Maximum 5 VDC
Ethernet port	IEEE P802.3


Electrical connections and specifications

Parameter	Description
Power source voltage	100–240 V ± 10%
Power frequency	50–60 Hz ± 10%
Power used	Maximum 600 VA
Power cord length	2.5 m plug: IEC13 female
Operational temperature	4–40 °C
Altitude	Up to 3000 m
RFID frequency	13.56 MHz
Tank material	High purity polyethylene

Weights

		Dry weight	Shipping weight	Operating weight
Purification unit	7003/7005	21.0 kg	24.0 kg	26.0 kg
	7010	23.7 kg	26.9 kg	29.7 kg
	7015	24.0 kg	27.2 kg	30.0 kg
Storage tank	25 L	6.7 kg	8.5 kg	31.7 kg
	50 L	7.6 kg	10.6 kg	57.6 kg
	100 L	10.9 kg	12.8 kg	110.9 kg

E-POD® and Q-POD® dispenser specifications

Tubing and port requirements

Parameter	Description
Dispenser tubing length	0.92 m
Distance from purification unit to POD dispenser	Maximum 5 m
Distance between two POD dispensers (Maximum 4 dispensers connected in series)	Maximum 5 m
POD dispenser data connection with unit	Ethernet
Electrical connection	Powered by purification unit (24–28 VDC)
Foot pedal port	3.3 V


Weights

	Dry weight	Shipping weight	Operating weight
Q-POD®/E-POD®	4.7 kg	7.2 kg	5.5 kg


Screen description and functionalities

Parameter	Description
Capacitive touchscreen	Size: 5 in; Resolution: 800 x 480
USB port	USB 2.0 Highspeed standard
Speaker	Impedance: 8 Ω /max output power: 0.5 W
Display in 9 languages	Chinese/English/French/German/Italian/Japanese/ Portuguese/Russian/Spanish


Examples of installation options


Up to 4 PODs connected

Under-bench installation

Bench-top and under-bench installation

Wall-mounted installation

Water specifications international conformity

Feed water requirements

Feed water	Potable tap water
Pressure	1–6 bar
Temperature	5–35 °C
Conductivity	<2000 µS/cm at 25 °C
Dissolved CO ₂	<30 ppm
Free chlorine	<3 ppm
Fouling Index	<10
pH	4–10
Total Organic Carbon (TOC)	<1 ppm
Langlier Saturation Index (LSI)	<0.3
Hardness (as CaCO ₃)	<300 ppm
Silica	<30 ppm

Ultrapure, Type 1 water specifications¹ (from Q-POD® dispenser)

Resistivity at 25 °C ²	18.2 MΩ · cm
Conductivity at 25 °C	0.055 µS/cm
TOC	≤2 ppb (µg/L) ³ ; typically ≤5 ppb
Particles ⁴	No particles with size >0.22 µm
Bacteria ⁵	<0.01 CFU/mL
Pyrogens (endotoxins) ⁶	<0.001 EU/mL
RNases ⁷	<1 pg/mL
DNases ⁷	<5 pg/mL
Proteases ⁷	<0.15 µg/mL
Flow rate	Up to 2 L/min

1. These values are typical and may vary depending on the nature and concentration of contaminants in the feed water.
2. Resistivity can also be displayed non-temperature-compensated as required by USP.
3. In the appropriate operating conditions, otherwise typically ≤5 ppb.
4. With Millipak® or Millipak® Gold filter.
5. With Millipak® or Biopak® filter when installed and used in a laminar flow hood.
6. With Biopak® polisher when installed and used in a laminar flow hood.
7. With Biopak® polisher.

Pure, Type 2 water specifications¹

Resistivity at 25 °C ²	>5 MΩ · cm; typically 10–15 MΩ · cm
Conductivity at 25 °C	0.2 μS/cm; typically 0.1 μS/cm
TOC	≤30 ppb
Production flow rate	3 L/h (Milli-Q® IQ 7003) 5 L/h (Milli-Q® IQ 7005) 10 L/h (Milli-Q® IQ 7010) 15 L/h (Milli-Q® IQ 7015)

From an E-POD® dispenser with final filter, the following water quality specifications are achieved:¹

Particulates ³	No particles with size >0.22 μm
Bacteria ⁴	<0.01 CFU/mL
Pyrogens (endotoxins) ⁵	<0.001 EU/mL
RNases ⁶	<1 pg/mL
DNases ⁶	<5 pg/mL
Proteases ⁶	<0.15 μg/mL
Flow rate	Up to 2 L/min

1. These values are typical and may vary depending on the nature and concentration of contaminants in the feed water.
2. Resistivity can also be displayed non-temperature-compensated as required by USP.
3. With Millipak® or Millipak® Gold filter.
4. With Millipak® or Biopak® filter when installed and used in a laminar flow hood.
5. With Biopak® polisher when installed and used in a laminar flow hood.
6. With Biopak® polisher.

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정수 카트리지	카다로그 번호
Milli-Q® IQ 7003-5 purification kit	IQ700XPKT1
Milli-Q® IQ 7003-5 purification kit H (for hardness)	IQ700XPKTH1
Milli-Q® IQ 7010-15 purification kit	IQ70XXPKT1
Milli-Q® IQ 7010-15 purification kit H (for hardness)	IQ70XXPKTH1
Vent Filter HF (for high-flow applications)*	TANKVH1A1

Application POD-Paks	카다로그 번호
Millipak® 0.22 µm filter	MPGP002A1
Millipak® Gold 0.22 µm sterile filter	MPGPG02A1
Biopak® polisher	CDUFB10A1
LC-Pak® polisher	LCPAK00A1
EDS-Pak® polisher	EDSPAK0A1
VOC-Pak® polisher	V0CPAK0A1

부속물	카다로그 번호
System wall mounting bracket	SYSTFIXA1
Tank wall mounting bracket	TANKFIXA1
Foot pedal	ZMQSFTSA1
Water sensor	ZWATSENA1
External solenoid valve for feed water	EXTSV00A1
Connector 2 m system to POD	ZFC0NNS2Q
Connector 5 m system to POD	ZFC0NNS5Q
Connector 2 m system to storage tank	ZFC0NN2ST
Connector 5 m system to storage tank	ZFC0NN5ST

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This image shows a single sheet of white paper with horizontal ruling lines. The lines are evenly spaced and run across the width of the page. There are no margins, text, or other markings on the paper.

Milli-Q®

Lab Water Solutions

자세한 정보는,
[MerckMillipore.com/labwater](https://www.MerckMillipore.com/labwater)을 방문하십시오

© 2019 Merck KGaA, Darmstadt, Germany and/or its affiliates. All Rights Reserved. Merck, the vibrant M, Milli-Q, Q-POD, E-POD, Elix, ech₂o, IPAK Gard, IPAK Quanta, IQnano, A10, Millipak, BioPak, VOC-Pak, EDS-Pak and LC-Pak are trademarks of Merck KGaA, Darmstadt, Germany or its affiliates. All other trademarks are the property of their respective owners. Detailed information on trademarks is available via publicly accessible resources.

Lit. No. MK_BR2157KO