

# CULTURE AS NATURE INTENDED



PHO<sub>2</sub>X BOX

**BAKER**

*RUSKINN*

[www.bakerco.com](http://www.bakerco.com)

# PHO<sub>2</sub>X BOX - A NEW GAS CONTROLLED HYPOXIA SYSTEM

Baker Ruskinn's PhO<sub>2</sub>x Box is a new, easy to use and economical Physoxia/Hypoxia system designed in vitro cell culture experiments. PhO<sub>2</sub>x Box comprises a single Gas Controller (offering both O<sub>2</sub> and CO<sub>2</sub> control) and a Cell Culture Chamber. The Cell Culture Chamber can be placed on a lab bench, or be placed inside an Incubator or a Workstation.




Figure 1. Western blot from University of Oulu Finland, showing proof of hypoxic conditions in the PhO<sub>2</sub>x Box inside a standard CO<sub>2</sub> Incubator

\*Printed with permission



## Gas Controller

- Controls and Monitors O<sub>2</sub> and CO<sub>2</sub> levels with a single touchscreen
- Intuitive touchscreen with large font number display
- Small footprint (325mm wide x 298mm deep)
- Easy user set up, takes around 10 minutes
- Only needs N<sub>2</sub> and CO<sub>2</sub> cylinders for operation for lower running costs
- 2 year warranty for peace of mind and lower running costs

# PHO<sub>2</sub>X BOX INCLUDED FEATURES



## Touchscreen Control

- O<sub>2</sub> control (from 0.1% to 20.0% in 0.1% increments)
- CO<sub>2</sub> control (from 0.1% to 20.0% in 0.1% increments)
- Hypoxic Cycling
- Temperature Display (of Cell Culture Chamber)

## Data Log:

- Up to 12 Months data history
- One Data Set stored per minute, each Data Set comprises: Time, Date, O<sub>2</sub> (Set/Actual), CO<sub>2</sub> (Set/Actual)
- Stored on SD card provided

## Audible Alarms:

- Low Gas (either CO<sub>2</sub> or N<sub>2</sub>)

## Culture Chambers

Your choice of 4 types of Cell Culture Chamber. Each is gas tight, has removable shelving and is easy to clean

- Black (light reducing) - Small and Large
- Clear - Small and Large



## CULTURE CHAMBER DIMENSIONS

|                                     |  | SMALL<br>(Black or Clear)<br>Chamber                     |        | LARGE<br>(Black or Clear)<br>Chamber                        |        |
|-------------------------------------|--|--|--------|---|--------|
|                                     |  | mm   | inches | mm  | inches |
|                                     |  | External dimensions                                      | Width  | 355   | 14     |
| Height                              | 173                                    |  | 6.8    | 238   | 9.4    |
| Depth                               | 369                                    |  | 14.5   | 369   | 14.5   |
| Internal dimensions                 | Width                                  | 335  | 13.2   | 335   | 13.2   |
|                                     | Height                                 | 140  | 5.5    | 205   | 8.1    |
|                                     | Depth                                  | 280  | 11     | 280   | 11     |
|                                     | Volume                                 | 13.1 Litres  |        | 19.2 Litres   |        |
| Workstation capacity-96 well plates | Number of plates (128mm x 86mm x 17mm) | 8 plates on shelf, (4 stacks of 2 plates high on shelf)  |        | 24 plates, (4 stacks of 2 plates high per shelf, 3 shelves) |        |
| Workstation capacity-T 75 flasks    | Number of flasks (150mm x 80mm x 86mm) | 6 (on shelf)   |        | 18 (6 on each shelf)  |        |
| PhO <sub>2</sub> x Box System       | Kgs/Lbs                                | 5.5 Kg/12.1 lbs (chamber) & 5.5 kg/12.1 lbs (controller) |        | 6.5 kg/14.3 lbs (chamber) & 5.5 kg/12.1 lbs (controller)    |        |

Usable Internal Volume: 13.1 litres (Small Chamber) | 19.2 litres (Large Chamber)

| Gas Controller<br>External<br>dimensions |  | mm     | inches |      |
|--|--|--------|--------|------|
|  |  | Width  | 325    | 12.8 |
|  |  | Height | 238    | 9.4  |
|  |  | Depth  | 298    | 11.7 |

## PERFORMANCE DATA

|                        | Gas Type         |                   | Coefficient of Variation (based on >700 readings) |                   |
|------------------------|------------------|-------------------|---|-------------------|
|                        | O <sub>2</sub> % | CO <sub>2</sub> % | O <sub>2</sub> %                                  | CO <sub>2</sub> % |
| Settings on Controller | 0.1              | 5                 | <8%   | <5%               |
|                        | 1                | 5                 | <5%   | <5%               |
|                        | 5                | 5                 | <2%   | <5%               |

## PERFORMANCE DATA AT DIFFERENT GAS LEVEL SETTINGS

|     | 0.1% O <sub>2</sub> / 5.0% CO <sub>2</sub> |                 |
|-----|--|-----------------|
|     | O <sub>2</sub>                             | CO <sub>2</sub> |
|     | Mean                                       | 0.1             |
| %CV | 6.7  | 2.2             |

|     | 1% O <sub>2</sub> / 5.0% CO <sub>2</sub> |                 |
|-----|--|-----------------|
|     | O <sub>2</sub>                           | CO <sub>2</sub> |
|     | Mean                                     | 1.1             |
| %CV | 7.9                                      | 1.9             |

|     | 5% O <sub>2</sub> / 5.0% CO <sub>2</sub> |                 |
|-----|--|-----------------|
|     | O <sub>2</sub>                           | CO <sub>2</sub> |
|     | Mean                                     | 5.1             |
| %CV | 1.3                                      | 1.6             |

# OPTIONAL ACCESSORIES

Single Cable Gland - Must be factory fitted

Gas Sample Port

O<sub>2</sub> Meter

# POWER SPECIFICATIONS


Voltage: 100 Volts to 240 Volts AC

Power: 15 Watts


130kWh/Year (Based on 24 / 7 / 365)




Time to Set Point Large Culture Chamber O<sub>2</sub>-5% CO<sub>2</sub>-5%



30 Second Door Opening Recovery Large Culture Chamber O<sub>2</sub>-5% CO<sub>2</sub>-5%



5 Minute Door Opening Recovery Large Culture Chamber O<sub>2</sub>-5% CO<sub>2</sub>-5%





## CULTURE AS NATURE INTENDED

---

Our Baker Ruskinn products are packed with new, innovative features that allow you to study even the most complex cell interactions under perfect physiological oxygen conditions. Whether you're hoping to replicate the environment of blood vessels or lung tissue, Baker Ruskinn provides the best tools for the job.

**BAKER** RUSKINN

[www.bakerco.com](http://www.bakerco.com)

