
Operating Instructions

Biomedical Freezer

MDF-MU339

MDF-MU539

MDF-MU539

Please read the operating instructions carefully before using this product, and keep the operating instructions for future use.

See page 35 for model number.

CONTENTS

| | |
|--------------------------------------------------------------|----|
| INTRODUCTION | 3 |
| PRECAUTIONS FOR SAFE OPERATION | 4 |
| SYMBOLS ON UNIT | 9 |
| ENVIRONMENTAL CONDITIONS | 9 |
| FREEZER COMPONENTS | 10 |
| Control panel | 12 |
| INSTALLATION SITE | 13 |
| INSTALLATION | 14 |
| START-UP OF UNIT | 16 |
| TEMPERATURE SETTING | 17 |
| KEY LOCK FUNCTION | 17 |
| ALARM TEMPERATURE SETTING | 18 |
| SETTING OF ALARM DELAY TIME | 19 |
| SETTING OF ALARM RESUME TIME | 20 |
| SETTING OF ALARM BUZZER INTERLOCK MODE | 21 |
| CHANGE OF COMPRESSOR DELAY TIME | 22 |
| DURING/AFTER POWER FAILURE | 23 |
| REMOTE ALARM TERMINAL | 23 |
| ALARMS & SAFETY FUNCTIONS | 24 |
| ROUTINE MAINTENANCE | 25 |
| Cleaning the exterior, interior, and accessories | 25 |
| Defrosting | 25 |
| CALIBRATION | 26 |
| REPLACEMENT OF WORN-OUT PARTS | 26 |
| Replacement of the battery for the power-failure alarm | 26 |
| Replacement of the fan motor | 26 |
| TROUBLESHOOTING | 27 |
| DISPOSAL OF UNIT | 28 |
| Recycle of battery | 28 |
| Refrigerant | 28 |
| Decontamination of unit | 28 |
| TEMPERATURE RECORDER (OPTION) | 33 |
| SPECIFICATIONS | 34 |
| PERFORMANCE | 35 |
| SAFETY CHECK SHEET | 36 |

INTRODUCTION

- Read the operating instructions carefully before using the product and follow the instructions for safe operation.
- PHC Corporation takes no responsibility for safety if the product is not used as intended or is used with any procedures other than those given in the operating instructions.
- Keep the operating instructions in a suitable place so that they can be referred to as necessary.
- The operating instructions are subject to change without notice for improvement of performance or function.
- Contact our sales representative or agent if any page of the operating instructions is lost or the page order is incorrect, or if the instructions are unclear or inaccurate.
- No part of the operating instructions may be reproduced in any form without the express written permission of PHC Corporation.

IMPORTANT NOTICE

PHC Corporation. guarantees this product under certain warranty conditions. However, please note that PHC Corporation. Ltd. shall not be responsible for any loss or damage to the contents of the product.

<Intended Use>

This equipment is designed for low temperature storage of biomedical samples.

PRECAUTIONS FOR SAFE OPERATION

It is imperative that the user complies with the operating instructions as they contain important safety advice.

Items and procedures are described so that you can use this unit correctly and safely. Following these precautions will prevent possible injury to the user and any other person.

Precautions are illustrated in the following way:

WARNING

Warning indicates a potentially hazardous situation which, if not avoided, could result in serious injury or death.

CAUTION

Failure to observe CAUTION signs could result in injury to personnel and damage to the unit and associated property.

Symbols have the following meanings:

 This symbol means caution.

 This symbol means an action is prohibited.

 This symbol means an instruction must be followed.

Be sure to keep the operating instructions in a place that is accessible to users of this unit.

PRECAUTIONS FOR SAFE OPERATION

WARNING

-  **Do not use the unit outdoors.** Exposure to rain may cause leakage and/or electric shock.
-  **Only qualified engineers or service personnel should install the unit.** Installation by unqualified personnel may cause electric shock or fire.
-  **Install the unit in a location capable of bearing the total combined weight (product + optional accessories + stored items). After installing the unit, be absolutely sure to take precautions to prevent the unit from falling over.** If the unit is installed in a location which is not strong enough or if the proper precautions are not taken, the unit may fall over and cause injuries.
-  **Turn the leveling feet to separate the casters from the floor and secure the unit.** If they are left touching the floor, the unit may inadvertently move out of position when its door is opened or closed. It may cause injury.
-  **Do not install the unit where there are high levels of moisture or where it may be splashed with water.** Installing the unit where there are high levels of moisture or where it may be splashed with water may cause the insulation to deteriorate and give rise to leakage and/or electric shock.
-  **Do not install the unit in a location where flammable or volatile substances are present.** Installing the unit in a location where flammable or volatile substances are present may cause explosions and/or a fire.
-  **Do not install the unit in a location where corrosive gases such as acids are present.** Installing the unit in a location where corrosive substances are present may cause electrical components to corrode, leading to leakage and/or electric shock due to the deterioration of insulation resulting from corroded electrical components.
-  **Do not place this unit in a location where it is difficult to disconnect the power supply plug.** Failure to disconnect the power supply plug may cause fire in the event of a problem or malfunction.
-  **Be absolutely sure to earth (ground) the unit to prevent electric shock.** Failure to earth the product may give rise to electric shock. If necessary, ask a qualified contractor to do this work.
-  **Do not connect the earth wire to a gas pipe, water pipe or lightning rod when earthing the unit.** Earthing the unit improperly may give rise to electric shock.
-  **Connect the unit to a power source as indicated on the rating label attached to the unit.** Use of any other voltage or frequency other than that on the rating label may cause fire or electric shock.
-  **Never store volatile or flammable substances in this unit except in a sealed container.** Such substances may cause explosion or fire if they leak.
-  **Never insert metal objects such as pins and wires into any vent, gap, or outlet on the unit.** This may cause electric shock or injury by accidental contact with moving parts.

PRECAUTIONS FOR SAFE OPERATION

WARNING

-  **When handling harmful samples (for example, those which consist of toxic, pathogenic or radioactive substances), install the unit inside a designated isolation facility.** If the unit is installed in a location which is not an isolation facility, there may be detrimental effects on both people and the natural environment.
-  **Before proceeding with maintenance or checking the unit, and disconnect the power supply plug.** Performing the work while power is still flowing to the product or while the power supply plug is still connected may give rise to electric shock and/or injury.
-  **Do not touch any electrical parts (such as power supply plug) or operate switches with a wet hand.** This may cause electric shock.
-  **Wear protective gloves and mask during maintenance.** Touching or inhaling chemicals or aerosols from around the unit may be detrimental to health.
-  **Never splash water directly onto the unit** as this may cause electric shock or short circuit.
-  **Never put containers with liquid on top of the unit** as this may cause electric shock or short circuit if the liquid is spilled.
-  **Never damage the power supply cord or power supply plug (by breaking, adapting, placing near a source of heat, bending with force, twisting, pulling, adding weight, or binding).** A damaged power supply cord or power supply plug may cause electric shock, short circuit, or fire
-  **Never disassemble, repair, or modify the unit yourself.** A high-voltage area is located inside the unit. Any work carried out by unauthorized personnel may result in electric shock. Contact our sales representative or agent for maintenance or repair.
-  **Make sure the power supply plug is pushed fully in.** Faulty insertion of the power supply plug may cause electric shock or fire due to generation of heat. Never use a damaged power supply plug or loose power supply outlet
-  **Disconnect the power supply plug if there is anything wrong with the unit.** Continued abnormal operation may cause electric shock or fire.
-  **Grip the power supply plug when disconnecting the power supply cord from the outlet.** Pulling the power supply cord may cause electric shock or short circuit.
-  **Remove dust from the power supply plug periodically.** Dust on the power supply plug may cause insulation failure due to moisture and thus cause a fire. Disconnect the power supply plug and wipe it with a dry cloth

PRECAUTIONS FOR SAFE OPERATION

WARNING

-  **Disconnect the power supply plug before moving the unit.** Take care not to damage the power supply cord. A damaged power supply cord may cause electric shock or fire.
-  **Disconnect the power supply cord when the unit is not in use for long periods.** Keeping the unit connected may cause electric shock, leakage, or fire due to the deterioration of insulation.
-  If the unit is to be stored unused in an unsupervised area for a long period, **ensure that children do not have access and that doors cannot be closed completely.**
-  **Ask a qualified contractor to carry out disassembly and disposal of the unit.** Leaving the unit in a location that can be accessed by third parties may result in unexpected accidents (e.g. the unit may be used for unintended purposes).
-  **Do not leave the plastic bags used for packing in a place where they can be reached by small children** as this may result in unexpected accidents such as suffocation.
-  **Never replace the battery for the power-failure alarm yourself.** Only qualified engineers or service personnel should replace the battery.
-  **When moving the unit, be sure to take precautions to prevent it from falling over.** Moving the unit with too much force may cause it to fall over, possibly resulting in injury. A qualified individual must be assigned to supervise the safe movement and relocation of the unit.

PRECAUTIONS FOR SAFE OPERATION

CAUTION

-  **Never install the unit in a location where corrosive materials such as sulphur compounds are likely to be generated (e.g. near a drainage facility).** Corrosion of the copper pipes may result in the deterioration and consequently the failure of the cooling unit.
-  **This unit must be plugged into a dedicated circuit protected by a branch circuit-breaker.**
-  **Use a dedicated power source as indicated on the rating label attached to the unit.** A multiple-tap source may cause fire due to excessive heating.
-  **Do not climb on top of the unit or put any objects on the unit.** Falling from the unit may cause injury; falling objects may cause damage to the unit
-  **Never store corrosive substances such as acids or alkalis in this unit except in a sealed container.** These may be harmful to your health and may cause corrosion of internal components or electrical parts.
-  **Check the settings when restarting operation after a power failure or after turning the power off.** The settings may have changed as a result of stopping the unit. Stored items inside the unit may be adversely affected when operation is resumed if the settings have changed.
-  **To ensure the safety of the service engineer, submit a safety check sheet with the required items filled out.** This is provided as the photocopiable “Safety Check Sheet” at the end of these operating instructions
-  **Do not touch the power supply joint of the temperature recorder(option).** It may cause electric shock.
-  **To prevent frostbite, wear protective gloves when handling frozen items in the chamber.** Too much frost may cause chamber temperature rise resulting from incomplete door close.

SYMBOLS ON UNIT

The following symbols are attached to the unit. The table describes the meaning of the symbols.




| | |
|-----------------------------------------------------------------------------------|----------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------|
|  | This symbol is attached to covers that access high-voltage electrical components to prevent electric shock. Only a qualified engineer or service personnel should be allowed to open these covers. |
|  | This symbol indicates that caution is required. Refer to product documentation for details. |
|  | This symbol indicates an earth. |

ENVIRONMENTAL CONDITIONS

This equipment is designed to be safe at least under the following conditions (based on the IEC 61010-1):

- Indoor use;
- Altitude up to 2000 m;
- Ambient temperature 5 °C to 40 °C;
- Maximum relative humidity 80 % for temperature up to 31 °C decreasing linearly to 50 % relative humidity at 40 °C;
- Mains supply voltage fluctuations up to ± 10 % of the nominal voltage;
- Transient overvoltages up to the levels of OVERVOLTAGE CATEGORY II;
- Temporary OVERVOLTAGES occurring on the mains supply;
- Applicable pollution degree of the intended environment (POLLUTION DEGREE 2 in most cases)

FREEZER COMPONENTS


FREEZER COMPONENTS

- 1. Lock:** Turn counterclockwise to 180 ° with a key and the outer door is securely locked.
- 2. Door:** To open the door, grip the handle.
- 3. Handle:** Always grip this handle to open and close the outer door.
- 4. Space for temperature recorder:** A temperature recorder (optional component) can be attached here. See page 33 “TEMPERATURE RECORDER”.
- 5. Levelling feet <bottom>:** These are screw bolts used to install and fix the unit. Adjust the height of the levelling feet by turning the screw bolts until the 2 front casters are away from the floor. See page 14.
- 6. Caster:** 4 casters are provided to facilitate moving of the cabinet. For the installation, adjust the leveling feet so that the front 2 casters cannot contact with the floor.
- 7. Door latch:** To lock the outer door, turn this latch downward. To unlock, rotate this latch up. A padlock is also available.
- 8. Control panel:** To display the temperature setting and running condition. See page 12.
- 9. Drain hose <Right side>:** The drain hose can be used to collect the defrosted water. See page 25.
- 10. Fixture (on back side):** 2 fixtures are provided as spacers between the cabinet and wall and also serve as hooks to fix the unit. See page 14.
- 11. Power Supply cord:** When you connect the power supply plug to the power supply outlet, the main unit turns on.
- 12. Access port:** This is used for leading the measuring cable from the freezing chamber to the outside.
- 13. Battery switch:** Switch for battery used for power failure alarm. Always keep “ON”. Turn the switch “OFF” when the unit is in no use for a long period (more than 1 month).
- 14. Remote alarm terminal:** Used to notify an alarm condition of the unit to remote location. See page 23 “REMOTE ALARM TERMINAL”.

FREEZER COMPONENTS

Control panel

1. Alarm indicator (ALARM): This lamp is flashed when the audible alarm is activated.

2. Digital temperature indicator: Normally, the current chamber temperature is displayed in the temperature display. When the high/low temperature alarm is activated, the current chamber temperature blinks. When the self-diagnostic function detects an abnormality, the error code is displayed. See page 24.

3. Numerical value shift key (▲): Pressing this key in the setting mode causes the numerical value to shift. ON-OFF of key lock can be selected by pressing this key in the key lock mode.

4. Digit shift key (▶): Pressing this key in the setting mode causes the changeable digit to shift. Key lock mode is led by pressing this key for more than 5 seconds in the temperature display mode.

5. Set key (SET): Temperature setting mode is led by pressing this key. Once the key is pressed, the changeable digit is flashed. Pressing this key again after setting desired temperature, the setting is stored into computer memory. If there is no key operation for 90 seconds during the setting mode, the setting mode is invalid automatically.

6. Defrost key (DEF): When removing the frost, press this key for 5 seconds. The freezer operation is stopped. After removing the frost, press this key again. See page 25.

Note: The freezing operation never resumes automatically after defrosting.

7. Buzzer stop key (BUZZER): To silence the audible alarm, press this key.

8. Alarm test key (ALARM TEST): Test key for alarm device. By pressing this key, the alarm lamp is flashed, remote alarm is activated and buzzer sounds. This means all alarm function operate correctly.

INSTALLATION SITE

This unit must be installed in a location which meets all the conditions described below.

- If the unit is installed in a location which does not meet the conditions, its specified performance may not be achieved or malfunctions and accidents may occur.

■ A location not exposed to direct sunlight

Avoid any location which is exposed to direct sunlight. Installing the unit in a location exposed to direct sunlight may reduce its cooling performance.

■ A well-ventilated (airy) location

In order to ensure ventilation, leave clearances of at least 10 cm around the unit (at the left, right, top and back). Blocking the ventilation may reduce the unit's cooling performance or cause malfunctions.

■ A location away from sources of heat

Avoid any location which is close to a major source of heat (such as a heater or boiler). Installing the unit near a major source of heat may reduce the unit's cooling performance.

■ A location with minimal changes in temperature

Avoid any location where the ambient temperature is subject to sudden changes. If the unit is installed in a location where the ambient temperature is subject to sudden changes, it will not be possible to achieve a stable cooling performance. The ambient temperature to use is +5 °C degrees to +35 °C degrees for MDF-MU339, +5 °C degrees to +30 °C degrees for MDF-MU539.

■ A flat surface where the floor is also capable of bearing the total combined weight (product + optional accessories + stored items)

Install the unit on a flat surface which is even and which is capable of bearing the total combined weight (product + optional accessories + stored items). If the unit is installed where the surface is uneven or where the unit will be inclined at an angle, the unit will be unstable, and accidents or injuries may occur and/or unnecessary vibration or noise may be generated.

■ A location with minimal humidity

Install the unit in a location where the relative humidity is less than 80 %R.H. Installing the unit in a very humid location may cause earth faults and/or electric shock.

■ A location free of flammable or corrosive gases

Avoid any location exposed to flammable or corrosive gases. Flammable or corrosive gases can cause explosions and/or a fire. Furthermore, corrosion of the electrical parts may cause the insulation to be reduced and result in earth faults and/or electric shock.

■ A location where nothing can fall onto the unit

Avoid locations where objects may fall onto the unit. Objects falling and hitting the unit may cause it to break down or fail.

INSTALLATION

When installing the unit, follow the steps below to secure the unit properly, and also be absolutely sure to earth the unit.

1. Preparations after unpacking

Remove all the tape used to secure the doors and interior parts, and leave the doors open for a short while for ventilation.

If any surfaces of the outer cabinet are dirty, wipe the surface using a cloth moistened with a diluted neutral dish-washing detergent.

- Using an undiluted solution of detergent may cause the unit's plastic areas to crack. Follow the directions on the detergent for details of dilution.
- After wiping the unit using the diluted detergent, be absolutely sure to wipe the surfaces with a cloth dipped in clean water to remove traces of the detergent. After this, be absolutely sure to wipe the surfaces with a dry cloth, allowing the surfaces of the outer cabinet to dry out completely, and then proceed with the installation.

Note:


Remove the cable tie that bands the power supply cord. Prolonged contact with the tie may cause corrosion of the cord coating.

2. Securing and levelling the unit using the levelling feet

Rotate the front levelling feet clockwise until the casters are raised 5 mm to 10 mm above the floor surface [Fig. 1].

In addition, rotate the levelling feet slightly clockwise or anticlockwise, and adjust them so that the unit is completely level.


- When the casters are raised from the floor surface, the unit will be secured. If they are left touching the floor, the unit may accidentally move when its door is opened or closed.

3. Fixing the unit

Two fixtures are attached to the rear of the frame. Fix the frame to the wall with these fixtures and rope or chain.

(Fig. 2)

4. Ground (earth)

The ground (earth) is for preventing the electric shock in the case of the electrical insulation is somehow degraded. Always ground the unit at the time of installation.

- If the power supply outlet is not grounded, it will be necessary to install a ground by qualified engineers.

⚠ WARNING

Use a power supply outlet with ground (earth) to prevent electric shock. If the power supply outlet is not grounded, it is necessary to install a ground by qualified engineers.

Never ground the unit through a gas pipe, water main, telephone line or lightning rod. Such grounding may cause electric shock in the case of an incomplete circuit.


INSTALLATION

5. Installation of shelves

The shelves are packaged all together in the chamber. Fix the shelves on the shelf stoppers securely before start-up of the unit.

- Incomplete installation may cause injury or damage.

The shelf stoppers are attached to the standard location. Change the location of shelf stopper as necessary.

START-UP OF UNIT

Use the following procedure to start trial operation or actual operation of the unit.

1. In the empty state, power supply plug the plug into a power supply outlet. Confirm that the temperature display on the control panel is blinking.

2. Turn on the battery switch.

<Important>

Always turn the battery switch on when the unit is operating to ensure that the power-failure alarm is working.

3. Set the desired chamber temperature. See page 17 for the temperature setting.

4. Check that the chamber temperature reaches the desired temperature.

5. Make sure that the alarm lamp flashes and the buzzer sounds by pressing the alarm test key (ALARM TEST). The remote alarm is also operated. E09 is displayed on the digital temperature indicator and buzzer sounds if the battery switch is off. Make sure to turn on the battery switch.

6. After confirming the above, you can put articles into the biomedical freezer chamber in a small batch to prevent the temperature rise.

Note:

- When starting the operation of the biomedical freezer for the first time, the alarm lamp (ALARM) lights. When the chamber temperature reaches around the set temperature, then the alarm lamp goes out (The remote alarm is activated).
- If the battery switch is turned on before turning on the power of the biomedical freezer, the temperature alarm is activated and the buzzer sounds and the remote alarm is also activated. Check that the battery switch is off before turning on the biomedical freezer.

TEMPERATURE SETTING

Table 1 shows the basic procedure for setting the chamber temperature. Perform key operations in the sequence indicated in the table. The example in the table is based on the assumption that the desired temperature is -25 °C.

Note: The chamber temperature is set to -30 °C at the factory.

Table 1. Basic operation sequence (Example: Chamber temperature -25 °C)

| | Description of operation | Key operated | Indication after operation |
|---|--------------------------------------------------------------------|--------------|--------------------------------------------------------------------------------|
| 1 | Connect the power supply plug. | — | The current chamber temperature is displayed. |
| 2 | Press set key. | SET | The second digit is flashed. |
| 3 | Set to -25 with the numerical value shift key and digit shift key. | | When pressed, the settable digit is shifted. |
| | | | When pressed, the figure of settable digit changes. |
| 4 | Press set key. | SET | Set temperature is memorized and the current chamber temperature is displayed. |

- The chamber temperature can be set within the range of -35 °C to -18 °C. For MDF-MU339, the ambient temperature is 35 °C and the guaranteed temperature is -30 °C at no load. In the case of MDF-MU539, the ambient temperature is 30 °C and the guaranteed temperature is -30 °C at no load.
- The setting mode returns to the temperature display mode automatically when 90 seconds has passed without any key operation. In this case, any setting before pressing set key (SET) is not memorized.

KEY LOCK FUNCTION

This unit is provided with the key lock function. When the key lock is ON, change of temperature setting through the key pad is not possible. The key lock is set to OFF at the factory.

| Display | Mode | Function |
|---------|-----------------|---------------------------------------|
| L 0 | Key lock is OFF | Enable to change temperature setting |
| L 1 | Key lock is ON | Disable to change temperature setting |

Table 2. Procedure for key lock setting (change from key lock OFF to key lock ON)

| | Description of operation | Key operated | Indication after operation |
|---|-------------------------------------------------------------|--------------|--------------------------------------------------------------------------|
| 1 | — | — | The current chamber temperature is displayed. |
| 2 | Press digit shift key for 5 seconds. | | The display changes to L 0 and the first digit blinks. |
| 3 | Press numerical value shift key and scroll the figure to 1. | | When pressed, the figure of settable digit changes. |
| 4 | Press set key. | SET | The key lock is set to ON. The current chamber temperature is displayed. |

- The setting mode returns to the temperature display mode automatically when 90 seconds has passed without any key operation. In this case, any setting before pressing set key (SET) is not memorized.

ALARM TEMPERATURE SETTING

This unit is provided with both high and low temperature alarms. The temperature at which the alarm is activated can be changed. The high temperature alarm can be set to a temperature of + 5 °C to + 15 °C. The low temperature alarm can be set to a temperature between -15 °C and -5 °C.

Note: The temperature alarm is set at ± 10 °C of the set temperature at the factory.

As an example, Table 3 shows the procedure to set the high temperature alarm so that the alarm can activate when the chamber temperature is 5 °C higher than the set temperature.

Table 4 shows the procedure to set the low temperature alarm so that the alarm can activate when the chamber temperature is 5 °C lower than the set temperature.

Table 3. Procedure for setting high temperature alarm

| | Description of operation | Key operated | Indication after operation |
|---|------------------------------------------------------------------------------|--------------|----------------------------------------------------------------------------------|
| 1 | — | — | The current chamber temperature is displayed. |
| 2 | Press numerical value shift key for about 5 seconds. | ▲ | The display changes to F00 and the first digit blinks. |
| 3 | Set the figure to F01 by using digit shift key and numerical value shift key | ▲ | When pressed, the figure of settable digit changes. |
| 4 | Press set key. | SET | The current setting is displayed and the first digit is flashed. |
| 5 | Set the figure to 005 by using digit shift key and numerical value shift key | ▶ | When pressed, the settable digit is shifted. |
| | | ▲ | When pressed, the figure of settable digit changes. |
| 6 | Press set key. | SET | Alarm temperature is memorized and the current chamber temperature is displayed. |

Table 4. Procedure for setting low temperature alarm





| | Description of operation | Key operated | Indication after operation |
|---|------------------------------------------------------------------------------|--------------|----------------------------------------------------------------------------------|
| 1 | — | — | The current chamber temperature is displayed. |
| 2 | Press numerical value shift key for about 5 seconds. | ▲ | The display changes to F00 and the first digit blinks. |
| 3 | Set the figure to F02 by using digit shift key and numerical value shift key | ▲ | When pressed, the figure of settable digit changes. |
| 4 | Press set key. | SET | The current setting is displayed and the first digit is flashed. |
| 5 | Set the figure to -05 by using digit shift key and numerical value shift key | ▶ | When pressed, the settable digit is shifted. |
| | | ▲ | When pressed, the figure of settable digit changes. |
| 6 | Press set key. | SET | Alarm temperature is memorized and the current chamber temperature is displayed. |

• The setting mode returns to the temperature display mode automatically when 90 seconds has passed without any key operation. In this case, any setting before pressing set key (SET) is not memorized.

SETTING OF ALARM DELAY TIME

The delay time of the buzzer and remote alarm for high and low temperature alarm can be set between 0 and 15 minutes. The procedure in table 5 shows the sequence to set the delay time to 10 minutes. The delay time is set to 15 minutes at the factory.

Table 5 Setting procedure for alarm delay time (change from 15 minutes to 10 minutes)

| | Description of operation | Key operated | Indication after operation |
|---|------------------------------------------------------------------------------|--------------|--------------------------------------------------------------------------------------------------------------------------------------------------------------------|
| 1 | | ---- | The current chamber temperature is displayed.  |
| 2 | Press numerical value shift key for 5 seconds. | ▲ | The display changes to F00 and the first digit blinks..  |
| 3 | Set the figure to F50 by using digit shift key and numerical value shift key | ▶ | When pressed, the settable digit is shifted. |
| | | ▲ | When pressed, the figure of settable digit changes.  |
| 4 | Press set key. | SET | The current delay time is displayed. The first digit is flashed.  |
| 5 | Set the figure to 010 with the numerical value shift key. | ▲ | When pressed, the figure of settable digit changes.  |
| 6 | Press set key. | SET | The delay time is memorized and the current chamber temperature is displayed.  |

- The alarm delay time set mode returns to the temperature display mode automatically when 90 seconds has passed without any key operation.

SETTING OF ALARM RESUME TIME

The buzzer is cancelled by pressing buzzer stop key (BUZZER) on the control panel during alarm condition. However, if the alarm condition is continued after the “alarm resume time” has passed, the buzzer sounds again.

It is possible to change the alarm resume time by the procedure shown in the Table.6.

The example in the table is based on the assumption that the desired duration is 20 minutes.

Note: Factory setting; 30 minutes.

Table 6. Setting procedure for alarm resume time (change from 30 minutes to 20 minutes)

| | Description of operation | Key operated | Indication after operation |
|---|------------------------------------------------------------------------------|--------------|----------------------------------------------------------------------------|
| 1 | — | — | The current chamber temperature is displayed. |
| 2 | Press digit shift key for 5 seconds. | ▲ | The display changes to F00 and the first digit blinks.. |
| 3 | Set the figure to F25 by using digit shift key and numerical value shift key | ▶ | When pressed, the settable digit is shifted. |
| | | ▲ | When pressed, the figure of settable digit changes. |
| 4 | Press set key. | SET | The current resume time is displayed. The second digit is flashed. |
| 5 | Set the figure to 020 with the numerical value shift key. | ▲ | When pressed, the figure of settable digit changes. |
| 6 | Press set key. | SET | The setting is memorized and the current chamber temperature is displayed. |

- The settable alarm resume time are 10, 20, 30, 40, 50, or 60 minutes. The buzzer would not resume if the resume time is set in 000.
- The setting mode returns to the temperature display mode automatically when 90 seconds has passed without any key operation. In this case, any setting before pressing set key (SET) is not memorized.

SETTING OF ALARM BUZZER INTERLOCK MODE

Setting the alarm buzzer interlock mode enables the alarm buzzer and the remote buzzer to operate in conjunction / non-interlocking. Setting the alarm buzzer interlock mode is as follows in Table 7.

Note: Factory setting; non-interlocking "000".

Table 7. Change procedure for Interlock / non-interlock (changed from "000" to "001")

| | Description of operation | Key operated | Indication after operation |
|---|------------------------------------------------------------------------------|--------------|-------------------------------------------------------------------------------|
| 1 | | ---- | The current chamber temperature is displayed. |
| 2 | Press numerical value shift key for 5 seconds. | ▲ | The display changes to F00 and the first digit blinks. |
| 3 | Set the figure to F24 by using digit shift key and numerical value shift key | ▶ | When pressed, the settable digit is shifted. |
| | | ▲ | When pressed, the figure of settable digit changes. |
| 4 | Press set key. | SET | The current delay time is displayed. The first digit is flashed. |
| 5 | Set the figure to 001 with the numerical value shift key. | ▲ | When pressed, the figure of settable digit changes. |
| 6 | Press set key. | SET | The delay time is memorized and the current chamber temperature is displayed. |

- The alarm delay time set mode returns to the temperature display mode automatically when 90 seconds has passed without any key operation.

CHANGE OF COMPRESSOR DELAY TIME

The delay time of the compressor can be changed to reduce the load on the power line and to facilitate the start-up (reset) of the freezer after power failure.

The example in the table is based on the assumption that the delay time is changed to 4 minutes. The delay time is set in 3 minutes at the factory.

Note:

- The setting range for delay time is between 3 and 15 minutes. The cool down of chamber temperature may be slow when the setting of delay time is over 5 minutes, depending on the installation environment. There is no need of changing the delay time when the capacity of power source is adequate.

Table 8. Changing procedure for delay time (change from 3 minutes to 4 minutes)

| | Description of operation | Key operated | Indication after operation |
|---|------------------------------------------------------------------------------|--------------|-------------------------------------------------------------------------------|
| 1 | | ---- | The current chamber temperature is displayed. |
| 2 | Press numerical value shift key for 5 seconds. | ▲ | The display changes to F00 and the first digit blinks. |
| 3 | Set the figure to F05 by using digit shift key and numerical value shift key | ▲ | When pressed, the figure of settable digit changes. |
| 4 | Press set key. | SET | The current delay time is displayed. The first digit is flashed. |
| 5 | Set the figure to 004 with the numerical value shift key. | ▲ | When pressed, the figure of settable digit changes.. |
| 6 | Press set key. | SET | The delay time is memorized and the current chamber temperature is displayed. |

- The compressor starts to run with the delay time when the power supply cord is connected to the outlet or after power failure.
- The compressor delay time set mode returns to the temperature display area mode automatically when 90 seconds has passed without any key operation.

DURING/AFTER POWER FAILURE

■ Display of chamber temperature during power failure

By depressing the buzzer stop key (BUZZER) during “power failure alarm”, the buzzer stop and the chamber temperature is displayed on the temperature display for 5 seconds.

✧ Check the chamber temperature as appropriate by depressing the buzzer stop key (BUZZER) during power failure.

<Important>

The battery for power failure alarm is a consumable part. Replace the battery about every 3 years. The alarm would not be activated at power failure if the battery is not replaced regularly. Contact our sales representative or agent for the replacement of battery.

■ Operation check after recovery from the power failure

The set value is memorized by nonvolatile memory. Accordingly, the freezer resumes the operation with setting before power failure.

When the freezer is recovered from power failure with the chamber temperature higher than the preset temperature, then the high temperature alarm is activated and the buzzer sounds and the remote alarm is also activated. Please push the buzzer stop key (BUZZER) to silence buzzer and take appropriate actions if needed.

REMOTE ALARM TERMINAL


The alarm is relayed to a remote location when a remote alarm device (commercially available) is connected to the remote alarm terminal. Installation of a remote alarm device is recommended when the unit is installed in an unattended location so that the operator of the unit is notified of the alarm.

■ Contact our sales representative or agent to arrange the installation of a remote alarm device.

The remote alarm terminal is installed at the back of the unit. The alarm is outputted from this terminal. Contact capacity is DC 30 V, 2 A.

Contact output:

| Terminal | Normal status | Abnormal status |
|-----------------------|---------------|-----------------|
| Between COM. and N.O. | open | close |
| Between COM. and N.C. | close | open |

Note:

- The remote alarm terminal will be in alarm mode when the power supply cord is unplugged because this will be interpreted as a power failure.
- It is recommended to use standard signal and interface cables with a maximum length of 30 meters.

ALARMS & SAFETY FUNCTIONS

This unit has the alarms and safety functions shown in Table 9, and also self diagnostic functions.

Table 9 Alarms and safety functions

| Alarm | Situation | Indication | Alarm buzzer | Safety functions |
|-----------------------------|---------------------------------------------------------------------------------------------------------------------------------------------------|----------------------------------------------------------------------------|------------------------------------------|-------------------------------------------------------------|
| High temperature alarm | If the chamber temperature is higher than the temperature at which the high temperature alarm is activated. | Alarm lamp is flashed. Digital temperature indicator is flashed. | Intermittent tone with 15 minutes delay. | Remote alarm with 15 minutes delay. |
| Low temperature alarm | If the chamber temperature is lower than the temperature at which the low temperature alarm is activated. | Alarm lamp is flashed. Digital temperature indicator is flashed. | Intermittent tone with 15 minutes delay. | Remote alarm with 15 minutes delay. |
| Power failure alarm | In the case of power failure. When the power to the unit is disconnected. | Alarm lamp is flashed. | Intermittent tone | Remote alarm. |
| Auto-return | When there is no key pressing in each setting mode for 90 seconds. | Chamber temperature is displayed. | — | Finishing of each Setting mode. |
| Key lock | When the key lock is ON. | — | — | Change of setting is disable. |
| Thermal sensor abnormality | If the thermal sensor is disconnected. | Alarm lamp is flashed. E01 and chamber temp. are displayed alternately. | Intermittent tone | Remote alarm. Continuous running. |
| | If the thermal sensor is short-circuited. | Alarm lamp is flashed. E02 and chamber temp. are displayed alternately. | Intermittent tone | Remote alarm. Continuous running. |
| Sensor abnormality | If the condenser sensor is disconnected. | Alarm lamp is flashed. E03 and chamber temp. are displayed alternately. | Intermittent tone | Remote alarm. |
| | If the condenser sensor is short-circuited. | Alarm lamp is flashed. E04 and chamber temp. are displayed alternately. | Intermittent tone | Remote alarm. |
| Battery switch check | When battery switch is OFF at the time of alarm test. | Alarm lamp is flashed. E09 is flashed. | Intermittent tone | Remote alarm. |
| Condenser temp. abnormality | In the event of failure of fan motor for cooling the compressor. When the ambient temperature exceeds the usable environmental condition, etc. | Alarm lamp is flashed. E10 and chamber temp. are displayed alternately. | Intermittent tone | Remote alarm. Compressor stop (At about 88 ° C. or more) |
| Battery check | When about 3 years has passed with the power ON. | F1 and chamber temp. are displayed alternately. | — | — |
| Fan motor check | When about 6 years has passed with the power ON. | F2 and chamber temp. are displayed alternately. | — | — |

Note:

- The above power failure alarm is available when the battery switch is on and the battery is charged. If the battery switch is off or the battery is discharged, only the remote alarm is activated.
- The power failure alarm can be kept about 12 hours with the battery charged fully. 2-day operation of the freezer is needed to charge the battery full.
- The chamber temperature is displayed for 5 seconds if the buzzer stop key (BUZZER) key is depressed during the power failure alarm. Then the buzzer is stopped. The alarm lamp keeps flashing.
- After power failure, the operation is resumed with the condition before power failure since the temperature setting and alarm temperature setting and alarm resume time are memorized in a nonvolatile memory.

ROUTINE MAINTENANCE

WARNING

Always disconnect the power supply to the unit prior to any repair or maintenance of the unit in order to prevent electric shock or injury.

Cleaning the exterior, interior, and accessories

Use a dry cloth to wipe down the outside and inside of the unit and all accessories. If the outside panels are dirty, clean them with a diluted neutral dish-washing detergent.

Wipe off condensation from the exterior of the cabinet with a dry, soft cloth.

- Using an undiluted solution of detergent may cause the unit's plastic areas to crack. Follow the directions on the detergent for details of dilution.
- After the wiping the cabinet or accessories with a diluted detergent, be absolutely sure to wipe the surfaces with a cloth dipped in clean water to remove traces of the detergent. After this, be absolutely sure to wipe the surfaces with a dry cloth.

<Important>

- Do not use a brush, an acid, a thinner, laundry soap, a powder detergent, or boiling water for cleaning. These may cause damage to painted surfaces or cause perishing of plastic and rubber components. Moreover, do not wipe plastic and rubber components with a volatile material.
- In order to maintain the unit's intended level of performance, always replace accessories that have been removed for cleaning.

Defrosting

This product is refrigerated by the direct cooling. When it uses for a long time, frost appears on the chamber wall. It cannot be cooled down when there is much amount of frost. The defrost is mentioned in the following.

Use the scraper provided for removing the frost if the biomedical freezer operation must be continued. Pay attention not to impact or damage the inner wall.

1. When defrosting, temporarily move all the contents of storage containers in the biomedical freezer to another low-temperature freezer.
2. Remove the storage container and the shelf.
3. Press defrost key (DEF) for 5 seconds to stop the freezing operation. While the freezing operation is stopped, the current chamber temperature and dF are displayed on the digital temperature display alternately.
4. After a several hours, check visually that all frost was removed completely.
5. Use the drain hose to drain accumulated water. When draining, remove the cap attached to the bottom of the cabinet and the hose
6. Wipe off moisture remaining in the cabinet and return the hose and drain cap to their original positions.
7. Restore the storage container and the shelf.
8. Press defrost key (DEF) so that the freezing operation can be started.
9. Once the chamber temperature has dropped to the desired temperature, place the original contents back in the biomedical freezer chamber.

Note:

After the defrosting, the freezing operation is never resumed automatically. Make sure to press defrost key (DEF) to start the freezing operation after defrosting.

CALIBRATION

During continuous operation, the following service tasks must be performed:

- Perform a temperature calibration at least once a year.

For temperature calibration, contact our sales representative or agent.

REPLACEMENT OF WORN-OUT PARTS

Replacement of the battery for the power failure alarm

Replace the battery for the power failure alarm about every 3 years to ensure the alarm will operate in the event of a power failure. Contact our sales representative or agent to arrange the replacement of the battery.

- The alarm function (blinking of alarm indicator, sounding of buzzer) will not operate when the battery for the power failure alarm is flat.
- The alarm indicator blinks and the buzzer sound as a result of power from the battery for the power failure alarm. The regular replacement of the battery for the power failure alarm is important to prevent an increase in chamber temperature in the event of unexpected circumstances.

<Important>

The used battery is a recyclable resource. Do not dispose of the battery. Always follow the procedure for recycling.

Replacement of the fan motor

Replace the fan motor for the cooling circuit about every 6 years to ensure the appropriate operation of the freezer. Contact our sales representative or agent to arrange the replacement of the fan motor for the cooling circuit.

- The appropriate operation of the unit is maintained by cooling the compressor with a fan. The motor to operate the fan (the fan motor for the cooling circuit) has a limited lifespan. A fan motor with diminished performance may lead to reduced cooling performance.

TROUBLESHOOTING

If the unit malfunctions, check out the following before calling for service.

<Attention>

If the malfunction is not resolved after checking the following items or if the malfunction is not shown in the table below, contact our sales representative or agent.

| Malfunction | Check/Remedy |
|-------------------------------------------------------------------|-----------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------|
| Nothing operates even when the power supply plug is plugged in | <ul style="list-style-type: none"> <input type="checkbox"/> The unit is not connected to the power supply properly. <input type="checkbox"/> The capacity and voltage of the power supply is not sufficient. <input type="checkbox"/> There is a power failure. <input type="checkbox"/> The circuit breaker on the supply circuit is activated. <input type="checkbox"/> The fuse on the supply circuit is blown. |
| The alarm is activated during operation | <ul style="list-style-type: none"> <input type="checkbox"/> The unit is not connected to the power supply properly. <input type="checkbox"/> The capacity and voltage of the power supply is not sufficient. <input type="checkbox"/> There is a power failure. <input type="checkbox"/> The circuit breaker on the supply circuit is activated. <input type="checkbox"/> The fuse on the supply circuit is blown. <input type="checkbox"/> The chamber temperature setting has been changed. <input type="checkbox"/> The door has been kept open for a long time. <input type="checkbox"/> Containers with a high temperature (load) have been put in the chamber. <input type="checkbox"/> The door is open. |
| No key operation is available | <ul style="list-style-type: none"> <input type="checkbox"/> The key lock is set in ON (L 1). → Set the key lock in OFF (L 0). |
| During setting mode, the mode returns to temperature display mode | <ul style="list-style-type: none"> ■ The unit will return from setting mode to temperature display mode automatically after 90 seconds if no key is operated (auto-return function). |
| Excessive noise | <ul style="list-style-type: none"> <input type="checkbox"/> The floor is not stable. <input type="checkbox"/> The installation site is not level. <input type="checkbox"/> The freezer is tilted. <input type="checkbox"/> The cabinet is touching the surrounding wall. |
| The chamber does not get cold enough | <ul style="list-style-type: none"> <input type="checkbox"/> Warm material has been put in the chamber. <input type="checkbox"/> The door is frequently opened. <input type="checkbox"/> The chamber temperature setting is too high. <input type="checkbox"/> The unit is in direct sunlight. <input type="checkbox"/> The equipment is not installed in the appropriate place described in this instruction manual <input type="checkbox"/> The ventilation around the unit is blocked. <input type="checkbox"/> There is a nearby heat source. <input type="checkbox"/> The ambient temperature is too high. <input type="checkbox"/> There are too many items stored inside the chamber. <input type="checkbox"/> The access port is not covered. → The access port should be covered with insulation and rubber caps when not in use. <input type="checkbox"/> The door seal is damaged. → If it is damaged, contact our sales representative or agent for replacement. <input type="checkbox"/> A foreign substance is located between door gaskets. |

DISPOSAL OF UNIT

WARNING

If the unit is to be stored unused in an unsupervised area for an extended period **ensure that children do not have access and doors cannot be closed completely.**

The disposal of the unit should be accomplished by appropriate personnel. Always remove doors to prevent accidents such as suffocation.

Recycle of battery

Ni-MH

- Label indication is obliged to comply with Japanese battery regulation.

- Label indication is obliged to comply with Taiwanese battery regulation.

Refrigerant

REGULATION (EU) No 517/2014 on Certain Fluorinated Greenhouse Gases

This product contains fluorinated greenhouse gases as refrigerant.

| Refrigerant | GWP | Refrigerant (net) (kg) | Weight of CO ₂ equivalent (ton) |
|-------------|------|------------------------|--------------------------------------------|
| R-404A | 3922 | MDF-MU339 : 0.186 | MDF-MU339 : 0.73 |
| | | MDF-MU539 : 0.227 | MDF-MU539 : 0.89 |

Decontamination of unit

Before disposing a biomedical freezer with biohazardous danger, decontaminate the biomedical freezer to the extent possible by the user

DISPOSAL OF UNIT

Disposal of Old Equipment and Batteries

Only for European Union and countries with recycling systems

These symbols on the products, packaging, and/or accompanying documents mean that used electrical and electronic products and batteries must not be mixed with general household waste.

For proper treatment, recovery and recycling of old products and used batteries, please take them to applicable collection points in accordance with your national legislation.

By disposing of them correctly, you will help to save valuable resources and prevent any potential negative effects on human health and the environment.

For more information about collection and recycling, please contact your local municipality.

Penalties may be applicable for incorrect disposal of this waste, in accordance with national legislation.

Note for the battery symbol (bottom symbol):

This symbol might be used in combination with a chemical symbol. In this case it complies with the requirement set by the Directive for the chemical involved.

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TEMPERATURE RECORDER (OPTION)

The chamber temperature can be monitored and recorded by installing an optional temperature recorder. For the installation of a temperature recorder, an optional recorder fixing is necessary.

- For the attachment of a temperature recorder, contact our sales representative or agent.
- For the usage of the temperature recorder, also refer to an installation sheet enclosed with the temperature recorder.

Main specifications of temperature recorder

| | MTR-G85C | MTR-4015LH |
|-------------------------------|--------------------------------------|---------------------------|
| Recording range | Between -100 °C and +40 °C | Between -40 °C and +14 °C |
| Feed speed of recording paper | 1-day/turn, 7-day/turn, 32-day/turn | 31-day/batch |
| Recording paper | Circular type | Strip type |
| Power source | Supplied from the biomedical freezer | Dry cell |
| Recorder fixing(Optional) | MPR-S7 | MPR-S30 |

SPECIFICATIONS

| Product name | Biomedical Freezer MDF-MU339 | Biomedical Freezer MDF-MU539 |
|------------------------|------------------------------------------------------------------------------------------------------------------------------------------------------------------------|---------------------------------|
| External dimensions | W616 mm x D770 mm x H1802 mm | W793 mm x D770 mm x H1802 mm |
| Internal dimensions | W472 mm x D614 mm x H1262 mm | W649 mm x D614 mm x H1262 mm |
| Effective capacity | 369 L | 504 L |
| Exterior | Painted steel | |
| Interior | Painted steel | |
| Insulation | Rigid polyurethane foamed-in place | |
| Outer door | Painted steel | |
| Shelf | Made of hard steel wire + PE coating : 6 pcs Load; 30 kg/shelf | |
| Evaporator | Tube on sheet type (Liner combined use) | |
| Access port | Diameter 30 mm, 1 on back side | |
| Condenser | Wire and tube type | |
| Compressor | Hermetic type, 250 W | |
| Refrigerant | R-404A | |
| Temperature controller | Microcomputer control system | |
| Temperature display | Digital display (between -50 °C and +50 °C) | |
| Temperature sensor | Thermistor sensor | |
| Temperature alarm | Flash of digital temperature indicator and alarm lamp, Buzzer, (Remote alarm) | |
| Remote alarm contact | Allowable contact capacity: DC 30 V, 2 A *1 | |
| Accessories | Key 1 set, Defrost spatula 1 pc | |
| | Storage container : 6 pcs | Storage container : 12 pcs |
| Weight | 135 kg | 163 kg |
| Battery | For power failure alarm, Nickel-metal-hydride battery, DC 6 V, 1100 mAh, Automatic charge | |
| Optional component | Temperature recorder (MTR-G85C)+Recorder fixing (MPR-S7), Temperature recorder (MTR-4015LH)+ Recorder fixing (MPR-S30) Interface board (MTR-480, MTR-L03) *1, *2 | |

*1: It is recommended to use standard signal and interface cables with a maximum length of 30 meters.

*2: For the data acquisition system MTR-5000 user only. Contact our sales representative or agent for purchase.

Note:

- Design or specifications are subject to change without notice.
- Refer to the updated catalogue when ordering an optional component.

PERFORMANCE

| | |
|------------------------------|------------------------------------------------------------------|
| Product name | Biomedical Freezer MDF-MU339 |
| Model number | MDF-MU339-PE |
| Cooling performance | -30 °C (ambient temperature; 35 °C, no load) |
| Temperature control range | -30 °C to -20 °C |
| Rated voltage | AC 220 V/230 V/240 V |
| Rated frequency | 50 Hz |
| Rated power consumption | 180 W/185 W/195 W |
| | Max. 260 W |
| Noise level | 42 dB [A] (background noise; 20 dB) |
| Maximum pressure | 2.6 MPa |
| Usable environment condition | Temperature; 5 °C to 35 °C Humidity; equal or less than 80 %R.H. |

- The unit with CE mark complies with EU directives.

| | |
|------------------------------|------------------------------------------------------------------|
| Product name | Biomedical Freezer MDF-MU539 |
| Model number | MDF-MU539-PE |
| Cooling performance | -30 °C (ambient temperature; 30 °C, no load) |
| Temperature control range | -30 °C to -20 °C |
| Rated voltage | AC 220 V/230 V/240 V |
| Rated frequency | 50 Hz |
| Rated power consumption | 190 W/195 W/205 W |
| | Max. 265 W |
| Noise level | 42 dB [A] (background noise; 20 dB) |
| Maximum pressure | 2.5 MPa |
| Usable environment condition | Temperature; 5 °C to 30 °C Humidity; equal or less than 80 %R.H. |

- The unit with CE mark complies with EU directives.

⚠ CAUTION

**Please fill in this form before servicing.
Hand over this form to the service engineer to keep for his and your safety.**

Safety check sheet

1. Freezer contents :

- Risk of infection: Yes No
Risk of toxicity: Yes No
Risk from radioactive sources: Yes No

(List all potentially hazardous materials that have been stored in this unit.)

Notes :

2. Contamination of the unit

Unit interior

- No contamination Yes No
Decontaminated Yes No
Contaminated Yes No

Others:

3. Instructions for safe repair/maintenance/disposal of the unit

- a) The unit is safe to work on Yes No
b) There is some danger (see below) Yes No

Procedure to be adhered to in order to reduce safety risk indicated in b) below.

Date :

Signature :

Address, Division :

Telephone :

| | | | |
|-------------------------------------|----------------|----------------|-----------------------|
| Product name: Biomedical Freezer | Model: MDF- | Serial number: | Date of installation: |
|-------------------------------------|----------------|----------------|-----------------------|

Please decontaminate the unit yourself before calling the service engineer.

MEMO

MEMO

MEMO

PHC Corporation

1-1-1 Sakada, Oizumi-machi, Ora-gun, Gunma 370-0596, Japan

© PHC Corporation 2018

Printed in Indonesia
LDCL045400-2
N0418-20619