# PHCbi

## Operating Instructions

Biomedical Freezer

## MDF-MU549DH

Please read the operating instructions carefully before using this product, and keep the operating instructions for future use.

See page 37 for model number.

# CONTENTS

INTRODUCTION .....	3
INTENDED USE AND PRECAUTIONS .....	3
EMC PERFORMANCE .....	3
PRECAUTIONS FOR SAFE OPERATION .....	4
Symbols on unit .....	4
LABELS ON UNIT .....	9
ENVIRONMENTAL CONDITIONS .....	9
FREEZER COMPONENTS	
Main body .....	10
Control panel .....	12
INSTALLATION SITE .....	13
INSTALLATION .....	14
START-UP OF UNIT .....	16
TEMPERATURE SETTING .....	17
KEY LOCK FUNCTION .....	18
ALARM TEMPERATURE SETTING .....	18
SETTING OF ALARM DELAY TIME .....	21
SETTING OF ALARM RESUME TIME .....	22
SETTING OF ALARM BUZZER LINK MODE .....	23
CHANGE OF COMPRESSOR DELAY TIME .....	24
DURING/AFTER POWER FAILURE .....	25
REMOTE ALARM TERMINAL .....	25
ALARMS & SAFETY FUNCTIONS .....	26
ROUTINE MAINTENANCE	
Cleaning the exterior, interior, and accessories .....	27
Defrosting .....	27
CALIBRATION .....	28
REPLACEMENT OF WORN-OUT PARTS	
Replacement of the battery for the power-failure alarm .....	28
TROUBLESHOOTING .....	29
DISPOSAL OF UNIT	
Recycle of battery .....	30
Decontamination of unit .....	30
TEMPERATURE RECORDER (OPTION) .....	35
STORAGE CONTAINER (OPTION) .....	35
SPECIFICATIONS .....	36
PERFORMANCE .....	37
SAFETY CHECK SHEET .....	38

# INTRODUCTION

- Read the operating instructions carefully before using the product and follow the instructions for safe operation.
- PHC Corporation takes no responsibility for safety if the product is not used as intended or is used with any procedures other than those given in the operating instructions.
- Keep the operating instructions in a suitable place so that they can be referred to as necessary.
- The operating instructions are subject to change without notice for improvement of performance or function.
- Contact our sales representative or agent if any page of the operating instructions is lost or the page order is incorrect, or if the instructions are unclear or inaccurate.
- No part of the operating instructions may be reproduced in any form without the express written permission of PHC Corporation.

## IMPORTANT NOTICE

PHC Corporation guarantees this product under certain warranty conditions. However, please note that PHC Corporation shall not be responsible for any loss or damage to the contents of the product.

# INTENDED USE AND PRECAUTIONS

This equipment is designed for low temperature storage of DNAs and plasma.

- The effective storage period depends on the sample condition and storage temperature. It is necessary to determine the storage temperature and period suitable for the purpose.
- For the plasma, the lower storage temperature should be required for longer term storage. It is recommended to store the plasma at -80 °C or lower.

# EMC PERFORMANCE

Emission: EN 61326-1  
Immunity: EN 61326-1

This product is intended for use in a basic electromagnetic environment.

# PRECAUTIONS FOR SAFE OPERATION

It is imperative that the user complies with the operating instructions as they contain important safety advice.

Items and procedures are described so that you can use this unit correctly and safely. Following these precautions will prevent possible injury to the user and any other person.

Precautions are illustrated in the following way:

## **WARNING**

Warning indicates a potentially hazardous situation which, if not avoided, could result in serious injury or death.

## **CAUTION**

Failure to observe CAUTION signs could result in injury to personnel and damage to the unit and associated property.

Symbols have the following meanings:

-  This symbol means caution.
-  This symbol means an action is prohibited.
-  This symbol means an instruction must be followed.

Be sure to keep the operating instructions in a place that is accessible to users of this unit.

## Symbols on unit

The following symbols are attached to the unit. The table describes the meaning of the symbols.

	This symbol is attached to covers that access high-voltage electrical components to prevent electric shock. Only a qualified engineer or service personnel should be allowed to open these covers.
	This symbol indicates that caution is required. Refer to product documentation for details.
	This symbol indicates an earth.
	This symbol indicates Incorrect usage could lead to a fire hazard.

# PRECAUTIONS FOR SAFE OPERATION

## **WARNING**

-  **Do not use the unit outdoors.** Exposure to rain may cause leakage and/or electric shock.
-  **Only qualified engineers or service personnel should install the unit.** Installation by unqualified personnel may cause electric shock or fire.
-  **Install the unit in a location capable of bearing the total combined weight (product + optional accessories + stored items). After installing the unit, be absolutely sure to take precautions to prevent the unit from falling over.** If the unit is installed in a location which is not strong enough or if the proper precautions are not taken, the unit may fall over and cause injuries.
-  **Turn the leveling feet to separate the casters from the floor and secure the unit.** If they are left touching the floor, the unit may inadvertently move out of position when its door is opened or closed. It may cause injury.
-  **Do not install the unit where there are high levels of moisture or where it may be splashed with water.** Installing the unit where there are high levels of moisture or where it may be splashed with water may cause the insulation to deteriorate and give rise to leakage and/or electric shock.
-  **Do not install the unit in a location where flammable or volatile substances are present.** Installing the unit in a location where flammable or volatile substances are present may cause explosions and/or a fire.
-  **Do not install the unit in a location where corrosive gases such as acids are present.** Installing the unit in a location where corrosive substances are present may cause electrical components to corrode, leading to leakage and/or electric shock due to the deterioration of insulation resulting from corroded electrical components.
-  **Do not place this unit in a location where it is difficult to disconnect the power supply plug.** Failure to disconnect the power supply plug may cause fire in the event of a problem or malfunction.
-  **Be absolutely sure to earth (ground) the unit to prevent electric shock.** Failure to earth the product may give rise to electric shock. If necessary, ask a qualified contractor to do this work.
-  **Do not connect the earth wire to a gas pipe, water pipe or lightning rod when earthing the unit.** Earthing the unit improperly may give rise to electric shock.
-  **Connect the unit to a power source as indicated on the rating label attached to the unit.** Use of any other voltage or frequency other than that on the rating label may cause fire or electric shock.
-  **Never store volatile or flammable substances in this unit except in a sealed container.** Such substances may cause explosion or fire if they leak.
-  **Never insert metal objects such as pins and wires into any vent, gap, or outlet on the unit.** This may cause electric shock or injury by accidental contact with moving parts.

# PRECAUTIONS FOR SAFE OPERATION

## **WARNING**

-  **When handling harmful samples (for example, those which consist of toxic, pathogenic or radioactive substances), install the unit inside a designated isolation facility.** If the unit is installed in a location which is not an isolation facility, there may be detrimental effects on both people and the natural environment.
-  **Before proceeding with maintenance or checking the unit, and disconnect the power supply plug.** Performing the work while power is still flowing to the product or while the power supply plug is still connected may give rise to electric shock and/or injury.
-  **Do not touch any electrical parts (such as power supply plug) or operate switches with a wet hand.** This may cause electric shock.
-  **Wear protective gloves and mask during maintenance.** Touching or inhaling chemicals or aerosols from around the unit may be detrimental to health.
-  **Never splash water directly onto the unit** as this may cause electric shock or short circuit.
-  **Never put containers with liquid on top of the unit** as this may cause electric shock or short circuit if the liquid is spilled.
-  **Never damage the power supply cord or power supply plug (by breaking, adapting, placing near a source of heat, bending with force, twisting, pulling, adding weight, or binding).** A damaged power supply cord or power supply plug may cause electric shock, short circuit, or fire.
-  **Never disassemble, repair, or modify the unit yourself.** A high-voltage area is located inside the unit. Any work carried out by unauthorized personnel may result in electric shock. Contact our sales representative or agent for maintenance or repair.
-  **Make sure the power supply plug is pushed fully in.** Faulty insertion of the power supply plug may cause electric shock or fire due to generation of heat. Never use a damaged power supply plug or loose power supply outlet.
-  **Disconnect the power supply plug if there is anything wrong with the unit.** Continued abnormal operation may cause electric shock or fire.
-  **Grip the power supply plug when disconnecting the power supply cord from the outlet.** Pulling the power supply cord may cause electric shock or short circuit.
-  **Do not pull the power supply cord, the power supply plug, or the product side plug (inlet part) hard.** If the cord breaks, it may cause electric shock or fire.

# PRECAUTIONS FOR SAFE OPERATION

## **WARNING**

-  **Remove dust from the power supply plug periodically.** Dust on the power supply plug may cause insulation failure due to moisture and thus cause a fire. Disconnect the power supply plug and wipe it with a dry cloth
-  **Disconnect the power supply plug before moving the unit.** Take care not to damage the power supply cord. A damaged power supply cord may cause electric shock or fire.
-  **Disconnect the power supply cord when the unit is not in use for long periods.** Keeping the unit connected may cause electric shock, leakage, or fire due to the deterioration of insulation.
-  **Do not use power cord for other electrical equipment.** Such power supply cord may cause fire or electric shock.
-  **Always use the removal power supply cord.** Other power supply cord may cause electric shock or fire.
-  If the unit is to be stored unused in an unsupervised area for a long period, **ensure that children do not have access and that doors cannot be closed completely.**
-  **Ask a qualified contractor to carry out disassembly and disposal of the unit.** Leaving the unit in a location that can be accessed by third parties may result in unexpected accidents (e.g. the unit may be used for unintended purposes).
-  **Do not leave the plastic bags used for packing in a place where they can be reached by small children** as this may result in unexpected accidents such as suffocation.
-  **Never replace the battery for the power-failure alarm yourself.** Only qualified engineers or service personnel should replace the battery.
-  **When moving the unit, be sure to take precautions to prevent it from falling over.** Moving the unit with too much force may cause it to fall over, possibly resulting in injury. A qualified individual must be assigned to supervise the safe movement and relocation of the unit.
-  **Flammable and explosive product.** The unit contains flammable refrigerant. When repairing or recycling, only trained service personnel will repair and follow the procedure below.
  - Well ventilate the room to prevent refrigerant accumulation.
  - Keep fire away when the refrigerant is contained in the product.
  - Do not damage or break the pipework.

# PRECAUTIONS FOR SAFE OPERATION

## CAUTION

-  **Never install the unit in a location where corrosive materials such as sulphur compounds are likely to be generated (e.g. near a drainage facility).** Corrosion of the copper pipes may result in the deterioration and consequently the failure of the cooling unit.
-  **This unit must be plugged into a dedicated circuit protected by a branch circuit-breaker.**
-  **Use a dedicated power source as indicated on the rating label attached to the unit.** A multiple-tap source may cause fire due to excessive heating.
-  **Do not climb on top of the unit or put any objects on the unit.** Falling from the unit may cause injury; falling objects may cause damage to the unit.
-  **Never store corrosive substances such as acids or alkalis in this unit except in a sealed container.** These may be harmful to your health and may cause corrosion of internal components or electrical parts.
-  **Check the settings when restarting operation after a power failure or after turning the power off.** The settings may have changed as a result of stopping the unit. Stored items inside the unit may be adversely affected when operation is resumed if the settings have changed.
-  **To ensure the safety of the service engineer, submit a safety check sheet with the required items filled out.** This is provided as the photocopyable “Safety Check Sheet” at the end of these operating instructions.
-  **Do not touch the power supply joint of the temperature recorder(option).** It may cause electric shock.
-  **To prevent frostbite, wear protective gloves when handling frozen items in the chamber.** Too much frost may cause chamber temperature rise resulting from incomplete door close.
-  **Use designated parts for parts replacement.** Using an incorrect part may cause fire.
-  **Do not give strong shock or vibration during movement or use.** The piping may be damaged, causing a fire.
-  **Flammable and explosive product.** The unit contains flammable refrigerant. Consult repair manual/owner's guide before attempting to install or service this product. All safety precautions must be followed.
-  **Flammable refrigerant used. Flammable and explosive product.** Dispose of properly in accordance with national regulations.

# LABELS ON UNIT

<Labels applied to the unit>

To avoid accidents, users are advised to read carefully the hazard labels found at key locations on the interior and exterior of the unit.

Possible Danger	Warning/Caution Type Location of Danger	Warning/Caution Label	Description of Danger
Personal injury Sample damage	<b>Frostbite</b> <b>Rise in chamber temperature</b> Interior		To prevent frostbite, wear protective gloves when handling frozen items in the chamber. Too much frost may cause chamber temperature rise resulting from incomplete door close.
Personal injury	<b>Electric shock</b> Electric box		Attached to covers that access high-voltage electrical components to prevent electric shock. Only a qualified engineer or service personnel should be allowed to open these covers.
Personal injury	<b>Flammable and explosive product</b> Interior		This product contains flammable refrigerant. Please follow the instructions when recycling.


# ENVIRONMENTAL CONDITIONS

This equipment is designed to be safe at least under the following conditions (based on the IEC 61010-1):


- Indoor use;
- Altitude up to 2000 m;
- Ambient temperature 5 °C to 40 °C;
- Maximum relative humidity 80 % for temperature up to 31 °C decreasing linearly to 50 % relative humidity at 40 °C;
- Mains supply voltage fluctuations up to  $\pm 10$  % of the nominal voltage;
- Transient overvoltages up to the levels of OVERVOLTAGE CATEGORY II;
- Temporary OVERVOLTAGES occurring on the mains supply;
- Applicable pollution degree of the intended environment (POLLUTION DEGREE 2 in most cases)

# FREEZER COMPONENTS


## Main body


<Front side>

<Back side>

<Right side>


# FREEZER COMPONENTS

- 1. Lock:** Turn counterclockwise to 180 ° with a key and the outer door is securely locked.
- 2. Door:** To open the door, grip the handle.
- 3. Handle:** Always grip this handle to open and close the door.
- 4-1. Space for temperature recorder in bottom compartment :** A temperature recorder (optional component) can be attached here. See page 35 “TEMPERATURE RECORDER”.
- 4-2. Space for temperature recorder in top compartment :** A temperature recorder (optional component) can be attached here. See page 35 “TEMPERATURE RECORDER”.
- 5. Levelling feet <bottom>:** These are screw bolts used to install and fix the unit. Adjust the height of the levelling feet by turning the screw bolts until the 2 front casters are away from the floor. See page 14 .
- 6. Caster:** 4 casters are provided to facilitate moving of the cabinet. For the installation, adjust the leveling feet so that the front 2 casters cannot contact with the floor.
- 7. Door latch:** To lock the outer door, turn this latch downward. To unlock, rotate this latch up. A padlock is also available.
- 8. Control panel:** To display the temperature setting and running condition. See page 12.
- 9. Fixture (on back side):** 2 fixtures are provided as spacers between the cabinet and wall and also serve as hooks to fix the unit. See page 14.
- 10. Access port:** This is used for leading the measuring cable from the freezing chamber to the outside.
- 11. Battery switch:** Switch for battery used for power failure alarm. Always keep “ON”. Turn the switch “OFF” when the unit is in no use for a long period (more than 1 month).
- 12. Remote alarm terminal:** Used to notify an alarm condition of the unit to remote location. See page 25 “REMOTE ALARM TERMINAL”.
- 13. Drain hose <Right side>:** The drain hose can be used to collect the defrosted water. See page 27.
- 14. Power supply cord (inlet type):** Do not disconnect the power supply plug or the product side plug (inlet part) carelessly. If the power supply cord is disconnected, the contents of the chamber may be deteriorated.

# FREEZER COMPONENTS

## Control panel

- 1. Alarm lamp (ALARM):** This lamp is flashed when the audible alarm is activated.
- 2. Bottom compartment lamp (BOTTOM):** This lamp is turned on when the bottom compartment is selected.
- 3. Top compartment lamp (TOP):** This lamp is turned on when the top compartment is selected.
- 4. Digital temperature indicator:** Normally, the current chamber temperature is displayed in the temperature display. When the high/low temperature alarm is activated, the current chamber temperature blinks. When the self-diagnostic function detects an abnormality, the error code is displayed. See page 26.
- 5. Numerical value shift key (▲):** Pressing this key in the setting mode causes the numerical value to shift. ON-OFF of key lock can be selected by pressing this key in the key lock mode.
- 6. Digit shift key (▶):** Pressing this key in the setting mode causes the changeable digit to shift. Key lock mode is led by pressing this key for more than 5 seconds in the temperature display mode.
- 7. Set key (SET):** Temperature setting mode is led by pressing this key. Once the key is pressed, the changeable digit is flashed. Pressing this key again after setting desired temperature, the setting is stored into computer memory. If there is no key operation for 90 seconds during the temperature setting mode, the temperature setting mode is invalid automatically.
- 8. Compartment select key (TOP/BOTTOM):** By pressing this key, internal temperature display of top compartment or bottom compartment and each setting is selected.
- 9. Defrost key (DEF):** When removing the frost, press this key for 5 seconds. The freezer operation is stopped. After removing the frost, press this key again. See page 27.  
**Note:** The freezing operation never resumes automatically after defrosting.
- 10. Alarm buzzer stop key (BUZZER):** To silence the audible alarm, press this key.
- 11. Alarm test key (ALARM TEST):** Test key for alarm device. By pressing this key, the alarm lamp is flashed, remote alarm is activated and buzzer sounds. This means all alarm function operate correctly.

# INSTALLATION SITE

This unit must be installed in a location which meets all the conditions described below.

- If the unit is installed in a location which does not meet the conditions, its specified performance may not be achieved or malfunctions and accidents may occur.

## ■ A location not exposed to direct sunlight

Avoid any location which is exposed to direct sunlight. Installing the unit in a location exposed to direct sunlight may reduce its cooling performance.

## ■ A well-ventilated (airy) location

In order to ensure ventilation, leave clearances of at least 10 cm around the unit (at the left, right, top and back). Blocking the ventilation may reduce the unit's cooling performance or cause malfunctions.

## ■ A location away from sources of heat

Avoid any location which is close to a major source of heat (such as a heater or boiler). Installing the unit near a major source of heat may reduce the unit's cooling performance.

## ■ A location with minimal changes in temperature

Avoid any location where the ambient temperature is subject to sudden changes. If the unit is installed in a location where the ambient temperature is subject to sudden changes, it will not be possible to achieve a stable cooling performance.

## ■ A flat surface where the floor is also capable of bearing the total combined weight (product + optional accessories + stored items)

Install the unit on a flat surface which is even and which is capable of bearing the total combined weight (product + optional accessories + stored items). If the unit is installed where the surface is uneven or where the unit will be inclined at an angle, the unit will be unstable, and accidents or injuries may occur and/or unnecessary vibration or noise may be generated.

## ■ A location with minimal humidity

Install the unit in a location where the relative humidity is less than 80 %R.H. Installing the unit in a very humid location may cause earth faults and/or electric shock.

## ■ A location free of flammable or corrosive gases

Avoid any location exposed to flammable or corrosive gases. Flammable or corrosive gases can cause explosions and/or a fire. Furthermore, corrosion of the electrical parts may cause the insulation to be reduced and result in earth faults and/or electric shock.

## ■ A location where nothing can fall onto the unit

Avoid locations where objects may fall onto the unit. Objects falling and hitting the unit may cause it to break down or fail.

# INSTALLATION

When installing the unit, follow the steps below to secure the unit properly, and also be absolutely sure to earth the unit.

## 1. Preparations after unpacking

Remove all the tape used to secure the doors and interior parts, and leave the doors open for a short while for ventilation.

If any surfaces of the outer cabinet are dirty, wipe the surface using a cloth moistened with a diluted neutral dish-washing detergent.

- Using an undiluted solution of detergent may cause the unit's plastic areas to crack. Follow the directions on the detergent for details of dilution.
- After wiping the unit using the diluted detergent, be absolutely sure to wipe the surfaces with a cloth dipped in clean water to remove traces of the detergent. After this, be absolutely sure to wipe the surfaces with a dry cloth, allowing the surfaces of the outer cabinet to dry out completely, and then proceed with the installation.

### Note:


Remove the cable tie that bands the power supply cord. Prolonged contact with the tie may cause corrosion of the cord coating.

## 2. Securing and levelling the unit using the levelling feet

Rotate the front levelling feet clockwise until the casters are raised 5 mm to 10 mm above the floor surface [Fig. 1].

In addition, rotate the levelling feet slightly clockwise or anticlockwise, and adjust them so that the unit is completely level.


- When the casters are raised from the floor surface, the unit will be secured. If they are left touching the floor, the unit may accidentally move when its door is opened or closed.

## 3. Fixing the unit

Two fixtures are attached to the rear of the unit. Fix the unit to the wall with these fixtures and rope or chain.

(Fig. 2)

## 4. Ground (earth)

The ground (earth) is for preventing the electric shock in the case of the electrical insulation is somehow degraded. Always ground the unit at the time of installation.

- If the power supply outlet is not grounded, it will be necessary to install a ground by qualified engineers.

### WARNING

**Use a power supply outlet with ground (earth)** to prevent electric shock. If the power supply outlet is not grounded, it is necessary to install a ground by qualified engineers.

**Never ground the unit through a gas pipe, water main, telephone line or lightning rod.** Such grounding may cause electric shock in the case of an incomplete circuit.


# INSTALLATION

## 5. Installation of shelves

The shelves are packaged all together in the chamber. Set the shelves firmly in place on the shelf stoppers at the standard locations. (Fig. 3)

### Note:

- Incomplete installation may cause injury or damage.

# START-UP OF UNIT

Use the following procedure to start trial operation or actual operation of the unit.

1. In the empty state, power supply plug the plug into a power supply outlet. Confirm that the temperature display on the control panel is blinking.

2. Turn on the battery switch.

**<Important>**

Always turn the battery switch on when the unit is operating to ensure that the power-failure alarm is working.

3. Set the desired chamber temperature. See page 17 for the temperature setting.

4. Check that the chamber temperature reaches the desired temperature.

5. Make sure that the alarm lamp flashes and the buzzer sounds by pressing the alarm test key (ALARM TEST). The remote alarm is also operated. E09 is displayed on the digital temperature indicator and buzzer sounds if the battery switch is off. Make sure to turn on the battery switch.

6. After confirming the above, you can put articles into the biomedical freezer chamber in a small batch to prevent the temperature rise.

**Note:**

■ When starting the operation of the biomedical freezer for the first time, the alarm lamp (ALARM) lights. When the chamber temperature reaches around the set temperature, then the alarm lamp goes out (The remote alarm is activated).

■ If the battery switch is turned on before turning on the power of the biomedical freezer, the temperature alarm is activated and the buzzer sounds and the remote alarm is also activated. Check that the battery switch is off before turning on the biomedical freezer.







# TEMPERATURE SETTING

Table 1 shows the basic procedure for setting the chamber temperature. Perform key operations in the sequence indicated in the table. The example in the table is based on the assumption that the desired temperature is -25 °C.

**Note:** The chamber temperature is set to -40 °C at the factory.

**Table 1. Basic operation sequence (Example: Chamber temperature -25 °C)**

	Description of operation	Key operated	Indication after operation
1	Connect the power supply plug.	-----	The current chamber temperature is displayed. 
2	By pressing the compartment select key, select the top compartment.	TOP BOTTOM	Top compartment lamp lights and the current chamber temperature of top compartment is displayed. 
3	Press set key.	SET	The second digit is flashed. 
4	Set to -25 with the numerical value shift key and digit shift key.	▶	When pressed, the settable digit is shifted.
		▲	When pressed, the figure of settable digit changes. 
5	Press set key.	SET	Set temperature is memorized and the current chamber temperature is displayed. 
6	By pressing the compartment select key, select the bottom compartment.	TOP BOTTOM	Top compartment lamp lights and the current chamber temperature of bottom compartment is displayed. 
7	Press set key.	SET	The second digit is flashed. 
8	Set to -25 with the numerical value shift key and digit shift key.	▶	When pressed, the settable digit is shifted.
		▲	When pressed, the figure of settable digit changes. 
9	Press set key.	SET	Set temperature is memorized and the current chamber temperature is displayed. 




- The chamber temperature can be set within the range of -45 °C to -18 °C, the ambient temperature is 35 °C and the guaranteed temperature is -40 °C at no load.
- The setting mode returns to the temperature display mode automatically when 90 seconds has passed without any key operation. In this case, any setting before pressing set key (SET) is not memorized.

# KEY LOCK FUNCTION

This unit is provided with the key lock function. When the key lock is ON, change of temperature setting through the key pad is not possible. The key lock is set to OFF at the factory.

Display	Mode	Function
L 0	Key lock is OFF	Enable to change temperature setting
L 1	Key lock is ON	Disable to change temperature setting

**Table 2. Procedure for key lock setting (change from key lock OFF to key lock ON)**

	Description of operation	Key operated	Indication after operation
1	-----	-----	The current chamber temperature is displayed. 
2	Press digit shift key for 5 seconds.	▶	The first digit is flashed. 
3	Set the figure to L 1 by using numerical value shift key.	▲	When pressed, the figure of settable digit changes. 
4	Press set key.	SET	The key lock is set to ON. The current chamber temperature is displayed. 

• The setting mode returns to the temperature display mode automatically when 90 seconds has passed without any key operation. In this case, any setting before pressing set key (SET) is not memorized.

# ALARM TEMPERATURE SETTING

This unit is provided with both high and low temperature alarms. The temperature at which the alarm is activated can be changed. The high temperature alarm can be set to a temperature of + 5 °C to + 15 °C. The low temperature alarm can be set to a temperature between -15 °C and -5 °C.

**Note:** The temperature alarm is set at  $\pm 10$  °C of the set temperature at the factory.

Display	Mode	Function
F01	High temperature alarm set for top	See Table 3 on page 19
F02	Low temperature alarm set for top	See Table 4 on page 19
F03	High temperature alarm set for bottom	See Table 5 on page 20
F04	Low temperature alarm set for bottom	See Table 6 on page 20

As an example, Table 3 and Table 4 shows the procedure to set the high temperature alarm so that the alarm can activate when the chamber temperature is 5°C higher than the set temperature.

Table 5 and Table 6 shows the procedure to set the low temperature alarm so that the alarm can activate when the chamber temperature is 5°C lower than the set temperature.

# ALARM TEMPERATURE SETTING

**Table 3 Procedure for setting high temperature alarm (Top compartment)**

	Description of operation	Key operated	Indication after operation
1	-----	-----	The current chamber temperature is displayed.
2	Press numerical value shift key for 5 seconds.	▲	The display changes to F00 and the first digit blinks.
3	Set the figure to F01 by using numerical value shift key.	▲	When pressed, the figure of settable digit changes.
4	Press set key.	SET	The current setting is displayed and the first digit is flashed.
5	Set the figure to 005 by using digit shift key and numerical value shift key.	▶	When pressed, the settable digit is shifted.
		▲	When pressed, the figure of settable digit changes.
6	Press set key.	SET	Alarm temperature is memorized and the current chamber temperature is displayed.

- The set mode returns to the temperature display mode automatically when 90 seconds has passed without any key operation. In this case, any setting before pressing set key (SET) is not memorized.

**Table 4 Procedure for setting low temperature alarm (Top compartment)**

	Description of operation	Key operated	Indication after operation
1	-----	-----	The current chamber temperature is displayed.
2	Press numerical value shift key for 5 seconds.	▲	The display changes to F00 and the first digit blinks.
3	Set the figure to F02 by using numerical value shift key.	▲	When pressed, the figure of settable digit changes.
4	Press set key.	SET	The current setting is displayed and the first digit is flashed.
5	Set the figure to -05 by using digit shift key and numerical value shift key.	▶	When pressed, the settable digit is shifted.
		▲	When pressed, the figure of settable digit changes.
6	Press set key.	SET	Alarm temperature is memorized and the current chamber temperature is displayed.

- The set mode returns to the temperature display mode automatically when 90 seconds has passed without any key operation. In this case, any setting before pressing set key (SET) is not memorized.

# ALARM TEMPERATURE SETTING

**Table 5 Procedure for setting high temperature alarm (Bottom compartment)**

	Description of operation	Key operated	Indication after operation
1	-----	-----	The current chamber temperature is displayed.
2	Press numerical value shift key for 5 seconds.	▲	The display changes to F00 and the first digit blinks.
3	Set the figure to F03 by using numerical value shift key.	▲	When pressed, the figure of settable digit changes.
4	Press set key.	SET	The current setting is displayed and the first digit is flashed.
5	Set the figure to 005 by using digit shift key and numerical value shift key.	▶	When pressed, the settable digit is shifted.
		▲	When pressed, the figure of settable digit changes.
6	Press set key.	SET	Alarm temperature is memorized and the current chamber temperature is displayed.

- The set mode returns to the temperature display mode automatically when 90 seconds has passed without any key operation. In this case, any setting before pressing set key (SET) is not memorized.

**Table 6 Procedure for setting low temperature alarm (Bottom compartment)**





	Description of operation	Key operated	Indication after operation
1	-----	-----	The current chamber temperature is displayed.
2	Press numerical value shift key for 5 seconds.	▲	The display changes to F00 and the first digit blinks.
3	Set the figure to F04 by using numerical value shift key.	▲	When pressed, the figure of settable digit changes.
4	Press set key.	SET	The current setting is displayed and the first digit is flashed.
5	Set the figure to -05 by using digit shift key and numerical value shift key.	▶	When pressed, the settable digit is shifted.
		▲	When pressed, the figure of settable digit changes.
6	Press set key.	SET	Alarm temperature is memorized and the current chamber temperature is displayed.

- The set mode returns to the temperature display mode automatically when 90 seconds has passed without any key operation. In this case, any setting before pressing set key (SET) is not memorized.

# SETTING OF ALARM DELAY TIME

The delay time of the buzzer and remote alarm for high and low temperature alarm can be set between 0 and 15 minutes. The procedure in table 7 shows the sequence to set the delay time to 10 minutes. The delay time is set to 15 minutes at the factory.

**Table 7. Setting procedure for alarm delay time (change from 15 minutes to 10 minutes)**

	Description of operation	Key operated	Indication after operation
1	-----	-----	The current chamber temperature is displayed. 
2	Press numerical value shift key for 5 seconds.	▲	The display changes to F00 and the first digit blinks. 
3	Set the figure to F50 by using digit shift key and numerical value shift key.	▶	When pressed, the settable digit is shifted.
		▲	When pressed, the figure of settable digit changes. 
4	Press set key.	SET	The current delay time is displayed. The first digit is flashed. 
5	Set the figure to 010 by using numerical value shift key.	▲	When pressed, the figure of settable digit changes. 
6	Press set key.	SET	The delay time is memorized and the current chamber temperature is displayed. 

- The alarm delay time set mode returns to the temperature display mode automatically when 90 seconds has passed without any key operation.

# SETTING OF ALARM RESUME TIME




The buzzer is cancelled by pressing buzzer stop key (BUZZER) on the control panel during alarm condition. However, if the alarm condition is continued after the “alarm resume time” has passed, the buzzer sounds again.

It is possible to change the alarm resume time by the procedure shown in the Table 8.

The example in the table is based on the assumption that the desired duration is 20 minutes.

**Note:** Factory setting; 30 minutes.

**Table 8. Setting procedure for alarm resume time (change from 30 minutes to 20 minutes)**

	Description of operation	Key operated	Indication after operation
1	-----	-----	The current chamber temperature is displayed. 
2	Press digit shift key for 5 seconds.	▲	The display changes to F00 and the first digit blinks. 
3	Set the figure to F25 by using digit shift key and numerical value shift key.	▶	When pressed, the settable digit is shifted.
		▲	When pressed, the figure of settable digit changes. 
4	Press set key.	SET	The current resume time is displayed. The second digit is flashed. 
5	Set the figure to 020 by using numerical value shift key.	▲	When pressed, the figure of settable digit changes. 
6	Press set key.	SET	The setting is memorized and the current chamber temperature is displayed. 





- The settable alarm resume time are 10, 20, 30, 40, 50, or 60 minutes. The buzzer would not resume if the resume time is set in 000.
- The setting mode returns to the temperature display mode automatically when 90 seconds has passed without any key operation. In this case, any setting before pressing set key (SET) is not memorized.

# SETTING OF ALARM BUZZER LINK MODE

Setting the alarm buzzer link mode enables the alarm buzzer and the remote buzzer to operate in conjunction / non- linking. Setting the alarm buzzer link mode is as follows in Table 9.

**Note:** Factory setting; non- linking “000”.

**Table 9. Change procedure for link / non- link (changed from “000” to “001”)**

	Description of operation	Key operated	Indication after operation
1	-----	-----	The current chamber temperature is displayed. 
2	Press numerical value shift key for 5 seconds.	▲	The display changes to F00 and the first digit blinks. 
3	Set the figure to F24 by using digit shift key and numerical value shift key.	▶	When pressed, the settable digit is shifted.
		▲	When pressed, the figure of settable digit changes. 
4	Press set key.	SET	The current setting is displayed and the first digit is flashed. 
5	Set the figure to 001 by using numerical value shift key.	▲	When pressed, the figure of settable digit changes. 
6	Press set key.	SET	The setting is memorized and the current chamber temperature is displayed. 

- The alarm buzzer link mode returns to the temperature display mode automatically when 90 seconds has passed without any key operation.

# CHANGE OF COMPRESSOR DELAY TIME




The delay time of the compressor can be changed to reduce the load on the power line and to facilitate the start-up (reset) of the freezer after power failure.

The example in the table is based on the assumption that the delay time is changed to 4 minutes. The delay time is set in 3 minutes at the factory.

**Note:**

- The setting range for delay time is between 3 and 15 minutes. The cool down of chamber temperature may be slow when the setting of delay time is over 5 minutes, depending on the installation environment. There is no need of changing the delay time when the capacity of power source is adequate.

**Table 10. Changing procedure for delay time (change from 3 minutes to 4 minutes)**

	Description of operation	Key operated	Indication after operation
1	-----	-----	The current chamber temperature is displayed. 
2	Press numerical value shift key for 5 seconds.	▲	The display changes to F00 and the first digit blinks. 
3	Set the figure to F05 by using numerical value shift key.	▲	When pressed, the figure of settable digit changes. 
4	Press set key.	SET	The current delay time is displayed. The first digit is flashed. 
5	Set the figure to 004 by using numerical value shift key.	▲	When pressed, the figure of settable digit changes. 
6	Press set key.	SET	The delay time is memorized and the current chamber temperature is displayed. 

- The compressor starts to run with the delay time when the power supply cord is connected to the outlet or after power failure.
- The compressor delay time set mode returns to the temperature display area mode automatically when 90 seconds has passed without any key operation.

# DURING/AFTER POWER FAILURE

## ■ Display of chamber temperature during power failure

By depressing the buzzer stop key (BUZZER) during “power failure alarm”, the buzzer stop and the chamber temperature is displayed on the temperature display for 5 seconds.

- Check the chamber temperature as appropriate by depressing the buzzer stop key (BUZZER) during power failure.

## <Important>

The battery for power failure alarm is a consumable part. Replace the battery about every 3 years. The alarm would not be activated at power failure if the battery is not replaced regularly. Contact our sales representative or agent for the replacement of battery.

## ■ Operation check after recovery from the power failure

The set value is memorized by nonvolatile memory. Accordingly, the freezer resumes the operation with setting before power failure.

When the freezer is recovered from power failure with the chamber temperature higher than the preset temperature, then the high temperature alarm is activated and the buzzer sounds and the remote alarm is also activated. Please push the buzzer stop key (BUZZER) to silence buzzer and take appropriate actions if needed.

# REMOTE ALARM TERMINAL


The alarm is relayed to a remote location when a remote alarm device (commercially available) is connected to the remote alarm terminal. Installation of a remote alarm device is recommended when the unit is installed in an unattended location so that the operator of the unit is notified of the alarm.

- Contact our sales representative or agent to arrange the installation of a remote alarm device.

The remote alarm terminal is installed at the back of the unit. The alarm is outputted from this terminal. Contact capacity is DC 30 V, 2 A.

Contact output:

Terminal	Normal status	Abnormal status
Between COM. and N.O.	open	close
Between COM. and N.C.	close	open

## Note:

- The remote alarm terminal will be in alarm mode when the power supply cord is unplugged because this will be interpreted as a power failure.
- It is recommended to use standard signal and interface cables with a maximum length of 30 meters.

# ALARMS & SAFETY FUNCTIONS

This unit has the alarms and safety functions shown in Table 11, and also self diagnostic functions.

**Table 11. Alarms and safety functions**

Indication		Situation	Alarm buzzer	Safety functions	Alarm	
Digital temp. indicator	Alarm lamp					
Digital temperature indicator is flashed.	Flashed	If the chamber temperature is higher than the temperature at which the high temperature alarm is activated.	Intermittent tone with 15 minutes delay.	Remote alarm with 15 minutes delay.	High temperature alarm	
		If the chamber temperature is lower than the temperature at which the low temperature alarm is activated.			Low temperature alarm	
-----		In the case of power failure. When the power to the unit is disconnected.	Intermittent tone	Remote alarm.	Power failure alarm	
Chamber temperature is displayed.	-----	When there is no key pressing in each setting mode for 90 seconds.	-----	Finishing of each Setting mode.	Auto-return	
-----		When the key lock is ON.		Change of setting is disable.	Key lock	
E01 and chamber temp. are displayed alternately.	Flashed	If the thermal sensor is disconnected.	Intermittent tone	Remote alarm. Continuous running.	Thermal sensor Abnormality (Top compartment)	
E02 and chamber temp. are displayed alternately.		If the thermal sensor is short-circuited.			Thermal sensor Abnormality (Bottom compartment)	
E03 and chamber temp. are displayed alternately.		If the thermal sensor is disconnected.		Remote alarm.	Sensor abnormality (Top compartment)	
E04 and chamber temp. are displayed alternately.		If the thermal sensor is short-circuited.			Sensor abnormality (Bottom compartment)	
E05 and chamber temp. are displayed alternately.		If the condenser sensor is disconnected.		Remote alarm.	Battery switch check	
E06 and chamber temp. are displayed alternately.		If the condenser sensor is short-circuited.			Condenser temp. Abnormality (Top compartment)	
E07 and chamber temp. are displayed alternately.		If the condenser sensor is disconnected.			Condenser temp. Abnormality (Bottom compartment)	
E08 and chamber temp. are displayed alternately.		If the condenser sensor is short-circuited.			Remote alarm. Compressor stop (At about 80 ° C. or more)	Fan motor Abnormality (Top compartment)
E09 is flashed.		When battery switch is OFF at the time of alarm test.				Fan motor Abnormality (Bottom compartment)
E10 and chamber temp. are displayed alternately.		In the event of failure of fan motor for cooling the compressor. When the ambient temperature exceeds the usable environmental condition, etc.		Remote alarm.	Circuit board Abnormality (Top compartment)	
E11 and chamber temp. are displayed alternately.		When an abnormality occurs in the rotational speed of the fan motor.			Circuit board Abnormality (Bottom compartment)	
E12 and chamber temp. are displayed alternately.					When an abnormality occurs in the PFC circuit board.	Circuit board Abnormality (Top compartment)
E13 and chamber temp. are displayed alternately.		Circuit board Abnormality (Bottom compartment)				
E14 and chamber temp. are displayed alternately.		When an abnormality occurs in the inverter circuit board.				Circuit board Abnormality (Top compartment)
E15 and chamber temp. are displayed alternately.				Circuit board Abnormality (Bottom compartment)		
E16 and chamber temp. are displayed alternately.						
E17 and chamber temp. are displayed alternately.						
F1 and chamber temp. are displayed alternately.	-----	When about 3 years has passed with the power ON.	-----	-----	Battery check	

Note:

- The above power failure alarm is available when the battery switch is on and the battery is charged. If the battery switch is off or the battery is discharged, only the remote alarm is activated.
- The power failure alarm can be kept about 12 hours with the battery charged fully. 2-day operation of the freezer is needed to charge the battery full.

# ROUTINE MAINTENANCE

## WARNING

Always disconnect the power supply to the unit prior to any repair or maintenance of the unit in order to prevent electric shock or injury.

## Cleaning the exterior, interior, and accessories

Use a dry cloth to wipe down the outside and inside of the unit and all accessories. If the outside panels are dirty, clean them with a diluted neutral dish-washing detergent.

Wipe off condensation from the exterior of the cabinet with a dry, soft cloth.

- Using an undiluted solution of detergent may cause the unit's plastic areas to crack. Follow the directions on the detergent for details of dilution.
- After the wiping the cabinet or accessories with a diluted detergent, be absolutely sure to wipe the surfaces with a cloth dipped in clean water to remove traces of the detergent. After this, be absolutely sure to wipe the surfaces with a dry cloth.

### <Important>

- Never pour water onto or into the unit. Doing so can damage the electric insulation and cause failure.
- Do not use a brush, an acid, a thinner, laundry soap, a powder detergent, or boiling water for cleaning. These may cause damage to painted surfaces or cause perishing of plastic and rubber components. Moreover, do not wipe plastic and rubber components with a volatile material.
- In order to maintain the unit's intended level of performance, always replace accessories that have been removed for cleaning.

## Defrosting

This product is refrigerated by the direct cooling. When it uses for a long time, frost appears on the chamber wall. It cannot be cooled down when there is much amount of frost. The defrost is mentioned in the following.

### ■ Defrosting by a defrost spatula

Use the defrost spatula provided for removing the frost if the biomedical freezer operation must be continued. Pay attention not to impact or damage the inner wall.

### ■ Stop operation and defrosting

Set top or bottom with the compartment select key and then do the following operation.

1. When defrosting, temporarily move all the contents of storage containers in the biomedical freezer to another low-temperature freezer.
2. Remove the storage container and the shelf.
3. Press defrost key (DEF) for 5 seconds to stop the freezing operation. While the freezing operation is stopped, the current chamber temperature and dF are displayed on the digital temperature display alternately.
4. After a several hours, check visually that all frost was removed completely.
5. Use the drain hose to drain accumulated water. When draining, remove the cap attached to the bottom of the cabinet and the hose
6. Wipe off moisture remaining in the cabinet and return the hose and drain cap to their original positions.
7. Restore the storage container and the shelf.
8. Press defrost key (DEF) so that the freezing operation can be started.
9. Once the chamber temperature has dropped to the desired temperature, place the original contents back in the biomedical freezer chamber.

### Note:

After the defrosting, the freezing operation is never resumed automatically. Make sure to press defrost key (DEF) to start the freezing operation after defrosting.

# CALIBRATION

During continuous operation, the following service tasks must be performed:

- Perform a temperature calibration at least once a year.

For temperature calibration, contact our sales representative or agent.

# REPLACEMENT OF WORN-OUT PARTS

## Replacement of the battery for the power failure alarm

Replace the battery for the power failure alarm about every 3 years to ensure the alarm will operate in the event of a power failure. Contact our sales representative or agent to arrange the replacement of the battery.

- The alarm function (blinking of alarm indicator, sounding of buzzer) will not operate when the battery for the power failure alarm is flat.
- The alarm indicator blinks and the buzzer sound as a result of power from the battery for the power failure alarm. The regular replacement of the battery for the power failure alarm is important to prevent an increase in chamber temperature in the event of unexpected circumstances.

### <Important>

The used battery is a recyclable resource. Do not dispose of the battery. Always follow the procedure for recycling.

# TROUBLESHOOTING

If the unit malfunctions, check out the following before calling for service.

## <Attention>

If the malfunction is not resolved after checking the following items or if the malfunction is not shown in the table below, contact our sales representative or agent.

Malfunction	Check/Remedy
Nothing operates even when the power supply plug is plugged in	<ul style="list-style-type: none"> <li><input type="checkbox"/> The unit is not connected to the power supply properly.</li> <li><input type="checkbox"/> The capacity and voltage of the power supply is not sufficient.</li> <li><input type="checkbox"/> There is a power failure.</li> <li><input type="checkbox"/> The circuit breaker on the supply circuit is activated.</li> <li><input type="checkbox"/> The fuse on the supply circuit is blown.</li> <li><input type="checkbox"/> The power supply plug (inlet part) on the product side is not connected properly.</li> </ul>
The alarm is activated during operation	<ul style="list-style-type: none"> <li><input type="checkbox"/> The unit is not connected to the power supply properly.</li> <li><input type="checkbox"/> The capacity and voltage of the power supply is not sufficient.</li> <li><input type="checkbox"/> There is a power failure.</li> <li><input type="checkbox"/> The circuit breaker on the supply circuit is activated.</li> <li><input type="checkbox"/> The fuse on the supply circuit is blown.</li> <li><input type="checkbox"/> The chamber temperature setting has been changed.</li> <li><input type="checkbox"/> The door has been kept open for a long time.</li> <li><input type="checkbox"/> Containers with a high temperature (load) have been put in the chamber.</li> <li><input type="checkbox"/> The door is open.</li> </ul>
No key operation is available	<ul style="list-style-type: none"> <li><input type="checkbox"/> The key lock is set in ON (L 1).</li> <li>→ Set the key lock in OFF (L 0).</li> </ul>
During setting mode, the mode returns to temperature display mode	<ul style="list-style-type: none"> <li>■ The unit will return from setting mode to temperature display mode automatically after 90 seconds if no key is operated (auto-return function).</li> </ul>
Excessive noise	<ul style="list-style-type: none"> <li><input type="checkbox"/> The floor is not stable.</li> <li><input type="checkbox"/> The installation site is not level.</li> <li><input type="checkbox"/> The freezer is tilted.</li> <li><input type="checkbox"/> The cabinet is touching the surrounding wall.</li> </ul>
The chamber does not get cold enough	<ul style="list-style-type: none"> <li><input type="checkbox"/> Warm material has been put in the chamber.</li> <li><input type="checkbox"/> The door is frequently opened.</li> <li><input type="checkbox"/> The chamber temperature setting is too high.</li> <li><input type="checkbox"/> The unit is in direct sunlight.</li> <li><input type="checkbox"/> The equipment is not installed in the appropriate place described in this instruction manual</li> <li><input type="checkbox"/> The ventilation around the unit is blocked.</li> <li><input type="checkbox"/> There is a nearby heat source.</li> <li><input type="checkbox"/> The ambient temperature is too high.</li> <li><input type="checkbox"/> There are too many items stored inside the chamber.</li> <li><input type="checkbox"/> The access port is not covered.</li> <li>→ The access port should be covered with insulation and rubber caps when not in use.</li> <li><input type="checkbox"/> The door seal is damaged.</li> <li>→ If it is damaged, contact our sales representative or agent for replacement.</li> <li><input type="checkbox"/> A foreign substance is located between door gaskets.</li> </ul>

# DISPOSAL OF UNIT

## WARNING

If the unit is to be stored unused in an unsupervised area for an extended period **ensure that children do not have access and doors cannot be closed completely.**

**The disposal of the unit should be accomplished by appropriate personnel. Always remove doors to prevent accidents such as suffocation.**

## Recycle of battery

Ni-MH

- Label indication is obliged to comply with Japanese battery regulation.

- Label indication is obliged to comply with Taiwanese battery regulation.

## Decontamination of unit

Before disposing a biomedical freezer with biohazardous danger, decontaminate the biomedical freezer to the extent possible by the user

# DISPOSAL OF UNIT

## Disposal of Old Equipment and Batteries

### Only for European Union and countries with recycling systems

These symbols on the products, packaging, and/or accompanying documents mean that used electrical and electronic products and batteries must not be mixed with general household waste.

For proper treatment, recovery and recycling of old products and used batteries, please take them to applicable collection points in accordance with your national legislation.

By disposing of them correctly, you will help to save valuable resources and prevent any potential negative effects on human health and the environment.

For more information about collection and recycling, please contact your local municipality.

Penalties may be applicable for incorrect disposal of this waste, in accordance with national legislation.

#### Note for the battery symbol (bottom symbol):

This symbol might be used in combination with a chemical symbol. In this case it complies with the requirement set by the Directive for the chemical involved.

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# TEMPERATURE RECORDER (OPTION)

The chamber temperature can be monitored and recorded by installing an optional temperature recorder. For the installation of a temperature recorder, an optional recorder fixing is necessary.

- For the attachment of a temperature recorder, contact our sales representative or agent.
- For the usage of the temperature recorder, also refer to an installation sheet enclosed with the temperature recorder.

Main specifications of temperature recorder

	MTR-G85C	MTR-4015LH
Recording range	Between -100 °C and +40 °C	Between -40 °C and +14 °C
Feed speed of recording paper	1-day/turn, 7-day/turn, 32-day/turn	31-day/batch
Recording paper	Circular type	Strip type
Power source	Supplied from the biomedical freezer	Dry cell
Recorder fixing *1 (For top side)	MDF-S740T	—
Recorder fixing *1 (For lower front)	MPR-S7	MPR-S30

\*1: Option

# STORAGE CONTAINER (OPTION)

	MDF-05SC
Inner dimensions	W280 mm x D552 mm x H157 mm
Material	ABS plastic
Package	2 pcs/set

- Be careful not to damage the storage container by the scraper when removing the frost.
- Be careful not to drop the storage container from the shelf during use.

# SPECIFICATIONS

Product name	Biomedical Freezer MDF-MU549DH
External dimensions	W793 mm x D770 mm x H1802 mm
Internal dimensions	W649 mm x D614 mm x H600 mm (2 chambers)
Effective capacity	479 L
Exterior	Painted steel
Interior	Painted steel
Insulation	Rigid polyurethane foamed-in place
Outer door	Painted steel
Shelf	Made of hard steel wire + PE coating : 6 pcs Load; 30 kg/shelf
Evaporator	Tube on sheet type (Liner combined use)
Access port	Diameter 30 mm, back side/ Each chamber
Condenser	Wire and tube type
Compressor	Hermetic type, 400 W/ Used in upper chamber and lower chamber respectively
Refrigerant	R-290
Temperature controller	Microcomputer control system
Temperature display	Digital display (between -50 °C and +50 °C)
Temperature sensor	Thermistor sensor/ Each chamber
Temperature alarm	Flash of digital temperature indicator and alarm lamp, Buzzer, (Remote alarm)
Remote alarm contact	Allowable contact capacity: DC 30 V, 2 A *1
Accessories	Key 1 set, Defrost spatula 1 pc
Weight	165 kg
Battery	For power failure alarm, Nickel-metal-hydride battery, DC 6 V, 1100 mAh, Automatic charge
Optional component	Temperature recorder (MTR-G85C)+Recorder fixing (MPR-S7),(MDF-S740T) Temperature recorder (MTR-4015LH)+Recorder fixing (MPR-S30) Interface board (MTR-480, MTR-L03) *1,*2 Storage container: 2 pcs/set (MDF-05SC)

\*1: It is recommended to use standard signal and interface cables with a maximum length of 30 meters.

\*2: For the data acquisition system MTR-5000 user only. Contact our sales representative or agent for purchase.

## Note:

- Design or specifications are subject to change without notice.
- Refer to the updated catalogue when ordering an optional component.

# PERFORMANCE

Product name	Biomedical Freezer MDF-MU549DH
Model number	MDF-MU549DH-PE
Cooling performance	-40 °C (ambient temperature; 35 °C, no load)
Temperature control range	-40 °C to -20 °C
Rated voltage	AC 220 V/230 V/240 V
Rated frequency	50 Hz
Rated power consumption	215 W
	Max.420 W
Noise level	42 dB [A] (background noise; 20 dB)
Maximum pressure	1.8 MPa
Usable environment condition	Temperature; 5 °C to 35 °C Humidity; equal or less than 80 %R.H.

- The unit with CE mark complies with EU directives.

**⚠ CAUTION**

**Please fill in this form before servicing.  
Hand over this form to the service engineer to keep for his and your safety.**

## Safety check sheet

1. Freezer contents :

- Risk of infection: Yes No  
Risk of toxicity: Yes No  
Risk from radioactive sources: Yes No

(List all potentially hazardous materials that have been stored in this unit.)

Notes :

2. Contamination of the unit

Unit interior

- No contamination Yes No  
Decontaminated Yes No  
Contaminated Yes No

Others:

3. Instructions for safe repair/maintenance/disposal of the unit

- a) The unit is safe to work on Yes No  
b) There is some danger (see below) Yes No

Procedure to be adhered to in order to reduce safety risk indicated in b) below.

Date :

Signature :

Address, Division :

Telephone :

Product name: Biomedical Freezer	Model: MDF-	Serial number:	Date of installation:
-------------------------------------	----------------	----------------	-----------------------

Please decontaminate the unit yourself before calling the service engineer.

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**PHC Europe B.V.**