

Avalo® Series Medical Cart Line

Critical Care Cart

The importance of efficient organization and plenty of storage capacity are self-evident when it comes to critical care carts.

Capsa Healthcare understands this and offers a full selection through our Critical Care Cart models. Choose from three standard models that feature an arrangement of 3", 6" and 10" drawers that can be custom divided to meet your precise storage needs. To ensure safety and security of your supplies, the Avalo Series Critical Care Cart features your choice of a proven and simple keyless entry system or Best® core-removable lock.



Accessories & Options



Bulk Storage Bins



Sharps Container




I.V. Pole



Waste Container



Tilt Bin Organizer



Slide-Out Shelf

Cart Specifications

	Surface Height	Depth	Width
Standard	43.5" / 111 cm	24" / 61 cm	31" / 79 cm
Intermediate	39.5" / 100 cm	24" / 61 cm	31" / 79 cm
Compact	36" / 91 cm	24" / 61 cm	31" / 79 cm

Drawer Specifications

	Height	Depth	Width
Three Inch	3" / 7.6 cm	17.5" / 45 cm	21.5" / 55 cm
Six Inch	6.25" / 16 cm	17.5" / 45 cm	21.5" / 55 cm
Ten Inch	10" / 25 cm	17.5" / 45 cm	21.5" / 55 cm

Model Options



Standard



Intermediate



Compact

Critical Care Cart Features

- Wide selection of accessories and organization systems
- Expansive work area with optional slide-out work surface
- Delivered fully assembled, no tools required
- Precision-welded steel frame with high-impact panels
- 5-year warranty
- Easy ordering, fast delivery

Specifications, availability & components are subject to change without notice.