

Avalo® Series Medical Cart Line Treatment Cart

The Avalo Series Treatment Cart offers huge capacity and flexible organization of treatment supplies in 3", 6", and 10" seamless drawers that feature a unique integrated divider system.

Stable performance and long-term value are assured by using high-density cart panels that prevent denting and corrosion. Every Capsa Healthcare Treatment Cart features an innovative keyless locking system that is as simple to use and maintain as a standard household flashlight.



Treatment Cart Features

- Wide selection of accessories and organization systems
- Expansive work area with optional slide-out work surface
- Delivered fully assembled, no tools required
- Precision-welded steel frame with high-impact panels
- 5-year warranty
- Easy ordering, fast delivery

Accessories & Options



Drawer Dividers



Sharps Container




I.V. Pole



Waste Container



Storage Module



Slide-Out Shelf

Cart Specifications

	Surface Height	Depth	Width
Standard	43.5" / 111 cm	24" / 61 cm	31" / 79 cm
Intermediate	39.5" / 100 cm	24" / 61 cm	31" / 79 cm
Compact	36" / 91 cm	24" / 61 cm	31" / 79 cm

Drawer Specifications

	Height	Depth	Width
Three Inch	3" / 7.6 cm	17.5" / 45 cm	21.5" / 55 cm
Six Inch	6.25" / 16 cm	17.5" / 45 cm	21.5" / 55 cm
Ten Inch	10" / 25 cm	17.5" / 45 cm	21.5" / 55 cm

Model Options



Standard



Intermediate



Compact

Specifications, availability & components are subject to change without notice.