




# M-Series Treatment Cart

# LTC Medication Handling Solutions

The Capsa **M3-TC Treatment Cart** offers lightweight performance, simple maneuverability, and large storage capacity. Smart design provides plenty of room for supplies with drawers that can be customized to ensure optimal organization. A perfect balance of form and function is assured with addition of a host of convenience accessories.



### Features & Options



Waste w/ Lid



Organizer Tray



5" Tracking Wheel



5" Braking Wheel



Writing Surface



Sharps Container



Chair Rail Bumpers



Supply Dividers



Drawer Bins

### Cart Specifications

**Cart Base** 39.4"(h) 26.7"(d) 29.4"(w)

### Drawer Specifications

Drawer	W x D (in.)	W x D (cm)	Height (in.)	Height (cm.)
3"	21.5 x 18.5	55 x 47	3.75	9.5
4"	22 x 18.5	56 x 47	4.25	10.8
7"	22 x 18.5	56 x 47	7	17.8
10"	22 x 18.5	56 x 47	10	25.4
11"	22 x 18.5	56 x 47	11.25	28.8

Specifications, availability & components are subject to change without notice.

### M3-TC Treatment Cart Features

- Integrated handle facilitates maneuverability
- Expansive work area with optional slide-out work surface
- Durable drawers with 3 different depths
- Full width drawer handle
- Full selection of convenience accessories & options
- Durable design to ensure long-term value
- Front footrest provides ergonomic support