

Force SpinPrep



The Force SpinPrep is a personal centrifuge that is ideal for drying microarray slides. The centrifuge was designed to be extremely easy to use. Simply close the lid and the rotor instantly reaches 4,800 rpm's. Press the lid release button to open the lid and the rotor quickly decelerates to a stop for removing samples.

- Designed to dry microarray slides in approximately 10 seconds
- Accommodates 2 standard size slides
- Compact design, small footprint
- Transparent lid for easy sample viewing

Ordering Information:

Cat. No.	Description
SBC150	Force SpinPrep with rotor and two slide cassettes, 115V
SBC150-2	Force SpinPrep with rotor and two slide cassettes, 230V EU
SBC150-3	Force SpinPrep with rotor and two slide cassettes, 230V UK
SBC150SC	Replacement slide cassettes, set of 2 for Force SpinPrep

