

SelectPette Multichannel Pipettes



- Volume lock mechanism
- Available in both eight and twelve channel models
- 360° rotating manifold
- Spring loaded shaft for easy tip fit
- Fully autoclavable

Specifications and Ordering Information:

Cat. No.	Volume Range	Channels	Accuracy	Precision	Tip
SBP108-10	1 - 10 μ L	8	\pm 8.0 to \pm 2.0	<6.0 to <1.2	10 μ L
SBP108-50	5 - 50 μ L	8	\pm 4.0 to \pm 1.6	<2.5 to <0.6	200 μ L
SBP108-200	20 - 200 μ L	8	\pm 3.0 to \pm 1.0	<1.5 to <0.6	200 μ L
SBP108-300	50 - 300 μ L	8	\pm 1.6 to \pm 1.0	<1.5 to <0.6	300 μ L
SBP112-10	1 - 10 μ L	12	\pm 8.0 to \pm 2.0	<6.0 to <1.2	10 μ L
SBP112-50	5 - 50 μ L	12	\pm 4.0 to \pm 1.6	<2.5 to <0.6	200 μ L
SBP112-200	20 - 200 μ L	12	\pm 3.0 to \pm 1.0	<1.5 to <0.6	200 μ L
SBP112-300	50 - 300 μ L	12	\pm 1.6 to \pm 1.0	<1.5 to <0.6	300 μ L

Accessories

Cat No.	Description
SBP989	Linear stand for six SelectPette variable volume single or multichannel pipettes
SBP986	Shelf clip for three SelectPette variable volume single or multichannel pipette

