

SelectPette Single Channel Pipettes



- Significantly reduced plunger forces
- Effortless volume adjustment with volume lock mechanism
- Superior accuracy and precision for critical applications
- Adjustable tip ejector for universal tip fit
- Fully autoclavable

Specifications and Ordering Information:

Cat. No.	Volume Range	Accuracy	Precision	Tip
SBP100-2	0.1 - 2µL	±12.0 to ±1.5	<6.0 to <0.7	10µL
SBP100-10	0.5 - 10µL	±4.0 to ±0.5	<2.8 to <0.4	10µL
SBP100-20	2 - 20µL	±3.0 to ±0.8	<1.5 to <0.3	200µL
SBP100-50	5 - 50µL	±2.5 to ±0.8	<2.0 to <0.4	200µL
SBP100-100	10 - 100µL	±1.6 to ±0.8	<0.8 to <0.2	200µL
SBP100-200	20 - 200µL	±1.2 to ±0.6	<0.6 to <0.2	200µL
SBP100-250	50 - 250µL	±1.0 to ±0.6	<0.4 to <0.3	200µL/300µL
SBP100-1000	100 - 1000µL	±1.6 to ±0.6	<0.4 to <0.2	1000µL
SBP100-5000	1000 - 5000µL	±0.6 to ±0.5	<0.3 to <0.2	5000µL
SBP100-10000	1000 - 10000µL	±2.5 to ±0.5	<0.6 to <0.2	10000µL

Accessories

Cat No.	Description
SBP985	Carousel stand for six SelectPette single channel pipettes
SBP989	Linear stand for six SelectPette variable volume single or multichannel pipettes
SBP986	Shelf clip for three SelectPette variable volume single or multichannel pipette

