


## Cryo Rack for Freezer



### 특징

- Stainless로 제작
- Slide Type으로 만들어 Rack을 다 꺼내지 않아도 Cryobox의 Sample을 꺼낼 수 있다.(SD type)

### Specifications

#### Upright Models

제품명	SD-212UK	SD-216UK	SD-220UK	SD-309UK	SD-312UK
MDF-U700VX	-	-	24	-	24
MDF-U500VX	-	-	16	-	16
MDF-U33V	-	12	-	-	12
MDF-U54V	-	-	16	-	16
MDF-U74V	-	-	24	-	24
MDF-U3386S	-	12	-	-	12
MDF-4186S	16	16	-	16	16
MDF-U5386S	-	-	16	-	16
MDF-U7386S	-	-	24	-	24

#### Chest Models

제품명	IR-207CK	IR-209CK	IR-210CK	IR-213CK	IR-305CK	IR-306CK	IR-307CK	IR-309CK
MDF-C8V1	6	-	-	-	6	-	-	-
MDF-193(AT)	6	-	-	-	6	-	-	-
MDF-594(AT)	-	-	-	24	-	-	-	24
MDF-794(AT)	-	-	-	36	-	-	-	36
MDF-1156(AT)	-	9	-	-	-	9	-	-
MDF-C2156(ATN)	-	-	15	-	-	-	15	-

### Ordering information

Model No.	Type	Box 용량	Size (WxDxH, mm)
SD-212UK	Upright	2" box 3 x 4	145 x 435 x 235
SD-216UK	Upright	2" box 4 x 4	145 x 580 x 235
SD-220UK	Upright	2" box 4 x 5	145 x 569 x 290
SD-309UK	Upright	3" box 3 x 3	145 x 435 x 246
SD-312UK	Upright	3" box 3 x 4	145 x 435 x 328
IR-207CK	Chest	2" box 7	144 x 144 x 400
IR-209CK	Chest	2" box 9	144 x 144 x 515
IR-210CK	Chest	2" box 10	144 x 144 x 572
IR-213CK	Chest	2" box 13	144 x 144 x 741
IR-305CK	Chest	3" box 5	144 x 144 x 410
IR-306CK	Chest	별도문의	별도문의
IR-307CK	Chest	별도문의	별도문의
IR-309CK	Chest	3" box 9	144 x 144 x 741