

MD Series Medication Cart

LTC Medication Dispensing Solutions



Artromick Series Medication Carts offer unparalleled performance in a timeless design to support a range of medication systems for long-term care settings. The Artromick MD Series offers the selection of 3 models to ensure the medication cart meets the capacity requirements of its application.

MD Series Features

Choose from 3", 6" or 9" drawers and/or a full range of cassette sizes will maximize your medication capacity and organization. Dual-sided patient drawer access supports the highest patient capacities and the most efficient drawer configurations.

- **Control Panel:** Standard on all models, the Control Panel organizes common supplies and provides additional shelving for bulk supplies or liquid storage.
- **Seamless Top Cap:** Expansive and durable, the seamless top cap is easy to clean and features an integrated handle and raised edge to contain spills.
- **Cart Height:** Complete flexibility assured with the choice of standard or compact height.
- **Cart Bumper:** All Artromick Medication Carts feature a non-marring full wraparound bumper.
- **Locking Bar:** Unique, pivoting locking bar means "one-step" sectional locking and the assurance of complete security.
- **Warranty:** Artromick's exclusive Medication Cart Protection Plan offers a 10 year limited warranty.

Accessories & Options




Drawer Dividers



Hard Top Surface




Multi-Cavity Set



Electronic Locking



Narcotic Box




Removable Trays



Writing Surface



Accessories Tray



Best Locks

Specifications (includes single end cabinet)

		Height	Depth	Width
S Model	Compact	38 7/8"	26 3/8"	26 1/8"
	Standard	41 7/8"	26 3/8"	26 1/8"
M Model	Compact	38 7/8"	26 3/8"	36 7/8"
	Standard	41 7/8"	26 3/8"	36 7/8"
L Model	Compact	38 7/8"	26 3/8"	44"
	Standard	41 7/8"	26 3/8"	44"

Specifications, availability & components are subject to change without notice.