



Bugbox Plus Specifications

External Dimensions	Width (side to side): 830 mm / 32.5"
	Depth: 660 mm / 26"
	Height (excluding stand): 650 mm / 25.5"
Internal Dimensions	Width: 500 mm / 19.5"
	Depth: 460 mm / 18"
	Height: 420 mm / 16.5"
Workstation Weight	approx: 55 kg / 121 lbs.
Interlock Dimensions	Width: 150 mm / 6"
	Depth: 230 mm / 9"
	Height: 190 mm / 7.5"
Interlock Cycle Time	35 seconds
Capacity (90-mm plates)	Maximum: 234
	Working: 180
	Interlock: 18