



Bugbox Specifications

External Dimensions	Width (side to side): 800 mm / 31.5"
	Depth: 660 mm / 26"
	Height (excluding stand): 650 mm / 25.5"
Internal Dimensions	Width: 500 mm / 19.5"
	Depth: 460 mm / 18"
	Height: 420 mm / 16.5"
Workstation Weight	approx: 45 kg / 99 lbs.
Interlock Dimensions	Width: 100 mm / 4"
	Depth: 100 mm / 4"
	Height: 200 mm / 8"
Interlock Cycle Time	15 seconds
Capacity (90-mm plates)	Maximum: 270
	Working: 200
	Interlock: 10