



CLEANING INSTRUCTIONS:



CLEAN THE JED™ WITH SOAP AND WATER OR AN ENZYMATIC CLEANER (E.G., EmPower®) TO REMOVE ANY STAIN OR DEBRIS. THEN DISINFECT WITH A HIGH-LEVEL DISINFECTANT (E.G., GERMICIDAL DISPOSABLE CLOTHS, Cidex®) AND ALLOW TO DRY BEFORE USE.

DO NOT SOAK IN ANY CLEANING SOLUTION. DO NOT PLACE ANY PART OF THE JED™ OR DISPOSABLE MANDIBLE CUPS IN AN AUTOCLAVE OR STERIS.