

Dry Bath Blocks

Construction	Anodized, high grade, non-porous aluminum
Dimensions (W x D)	
Single block	3.55 x 3.01 in/9 x 7.6 cm
Dual block	3.55 x 6.01 in/9 x 15.3 cm
Thermometer well	Available in all blocks except those used with plates



Ordering Information:

Cat. No.	Description
D1101*	Solid block for machining (no holes)
D1102*	Block, 48 x 0.2 mL PCR tubes or 6 x 0.2 mL strips
D1102A*	Block, 20 x 2.0 mL
D1105*	Block, 24 x 0.5 mL
D1105A*	Block, 24 x 1.5 mL
D1106*	Block, 35 x 6 mm
D1110*	Block, 20 x 10 mm
D1112*	Block, 20 x 12 mm
D1113*	Block, 20 x 13 mm
D1115-TALL*	Block, 12 x 15 mL centrifuge tubes
D1116*	Block, 12 x 15/16 mm
D1117*	Block, 12 x 17 mm
D1120*	Block, 6 x 20 mm
D1125*	Block, 6 x 25 mm
D1150-TALL*	Block, 5 x 50 mL centrifuge tubes
D1196-PCR	Single block, 96 well PCR plate, skirted or nonskirted (single block only)
D1296	Dual block, 96 well microtiter plate or 4 slides (dual block only)
D1296-PCR	Dual block, 96 well PCR plate, skirted or nonskirted (dual block only)
D12384	Dual block, 384 well PCR plate (dual block only)
D1100-PS	PopStopper for microtube blocks, D1102, D1102A, D1105 & D1105A
D1105A-RACK	AccuRack accessory for 20 x 1.5 mL tubes, compatible with D1105A

*Can be utilized with single, dual and four block units

