



SCI-tive Specifications

Workstation Fabrication	Solvent Bonded Acrylic
External Dimensions	Width (side to side, excludes Gas Mixer): 1660 mm / 65"
	Depth: 826 mm / 32.5"
	Height (excluding stand): 1077 mm / 42"
Internal Dimensions	Width: 1200 mm / 47"
	Depth: 600 mm / 23.5"
	Height: 690 mm / 27"
Workstation Weight	approx: 230 kg / 507 lbs.
Alarms	Temperature: Visual and Audible
	Humidity: Visual
	Gas Low Pressure: Visual and Audible
Interlock Dimensions	Width: 270 mm / 11"
	Depth: 200 mm / 8"
	Height: 210 mm / 8"
Interlock Cycle Time	60 seconds
Atmosphere Control	Temperature: Ambient +5°C - 45°C
	Humidity Control: Ambient-85% RH
	O2 Control: 0.0%-20.9%
	CO2 Control: 0.0%-30.0%